Chevalier Valvulotome

English - Instructions for Use

Chevalier Valvulotom

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Valvulotomo Chevalier

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Valvulotomo Chevalier

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Chevalier-valvulotoom

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Chevalier Valvulotome

Chevalier Valvulotome

(Model Numbers 40025, 40030, 40035)

English - Instructions for Use

STERILE	EO
----------------	-----------

Product Description

The Chevalier Valvulotome is a device used to disrupt the valves of veins which are to be used as grafts in the in situ or non-reversed bypass procedure. The Chevalier Valvulotome consists of a 70 cm stainless steel shaft with two tips:

- The cutting tip is shaped like an inverted tulip with four teeth at its base. This tip is designed to cut the vein valve cusps.
- The other tip is atraumatic and used to test the patency and the course of the vein before the valvulotomy.

The Chevalier Valvulotome is available in three tip diameters:

2.5 mm, 3.0 mm and 3.5 mm.

Indication for Use

The Chevalier Valvulotome is indicated for the treatment of peripheral arterial disease which includes popliteal artery, femoral artery, tibioperoneal trunk, posterior tibial artery, peroneal artery, and peripheral vascular disease which includes atherosclerosis.

Intended Purpose

The Chevalier Valvulotome is intended for the treatment of vascular disorders and more particularly for excising or disrupting venous valves.

Intended User

The Chevalier Valvulotome is a surgical tool intended for use by experienced vascular surgeons trained in the procedures for which they are intended.

Patient Population

Patients of any gender, age or ethnicity with vascular disorders. There is no data for the use of this device on pregnant women and children. It is at the surgeon's discretion on whether to use it on this population.

Part of the Body Contacted

The Chevalier Valvulotome will come in contact with veins.

Clinical Condition

In situ peripheral bypass, non-reversed translocated bypass.

Clinical Benefits

The clinical benefits associated with the use of the Chevalier Valvulotome are increased limb salvage, amputation-free survival, and overall survival, comparable to rates observed with similar Valvulotomes.

Contraindications

- Endarterectomy procedures
- Thrombolysis procedures
- Vein stripping procedures
- Embolectomy procedures
- Vessel dilation procedures

Warnings

1. The Chevalier Valvulotome should not be passed through the proximal anastomosis. If the proximal valves are close to the long saphenous vein arch disinsertion, they can be cut with scissors under visual control.
2. The Chevalier Valvulotome should be handled with care during the procedure to avoid any vein wall disruption or puncture.
3. Do not use if the package or device is damaged.
4. Do not use if inner packaging is opened outside a sterile environment.
5. The Chevalier Valvulotome is to be used by a qualified physician.
6. Do not pass the Valvulotome through a vessel that has undergone synthetic grafting or contains implants.
7. Do not use device for valvulotomy unless the target vein is fully distended by arterial blood flow (due to anastomosis of artery to vein) or by saline injection.

Precautions

- Do not use devices past the expiration date printed on the labeling.
- These devices are single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/Repackaging" section of this document for further details).
- Keep any dissected portion of the vein straight and avoid twisting the vein during valvulotomy.
- Use caution when tracking the Chevalier Valvulotome in thrombophlebitic veins.

Potential Complications

Potential Device-related Complications:

- Entrapment of cutting blades in branches or anastomoses
- Inadequate valve disruption
- Residual tributaries
- Vein Stricture

Potential Procedure-related Complications:

- | | | |
|-------------|---------------------------------|------------------------------|
| • Death | • Hemorrhage | • Thrombus formation |
| • Embolism | • Infections | • Toxicity |
| • Erythema | • Intimal disruption | • Vessel damage |
| • Hematoma | • Irritation | • Vessel Occlusion/ Stenosis |
| • Hemolysis | • Post-procedure wound necrosis | • Vessel wall perforation |

Opening The Package

The Chevalier valvulotome is packaged in double sealed pouches. Open the package observing the standard aseptic technique. Take the outer pouch in one hand and peel back. The inner pouch can be removed using sterile gloves or instruments.

How Supplied

The Chevalier Valvulotome is provided sterile and is intended for single-use only.

Instructions For Use

1. The Chevalier Valvulotome has been designed to perform valvulotomy of veins which are primed.
2. If necessary, the patency and the course of the vein should be tested with the atraumatic tip, before the valvulotomy.
3. After the proximal anastomosis has been performed, the correct size cutting tip is inserted distally and advanced proximally. It should not be passed through the proximal anastomosis.
4. The valvulotomy is performed by withdrawing the valvulotome back through the filled valve where a small resistance may be felt. The instrument is then passed back and forth through the valve a few more times, usually four or five, rotating a quarter of a turn each time until the valve cusps are destroyed.
5. If strong resistance is encountered, the valvulotome should be re-advanced, rotated a quarter of a turn and the pass attempt repeated.
6. The valves are deemed to be destroyed when pulsatile blood flow is observed distally. The distal anastomosis may then be performed.
7. Operative angiography is recommended for arteriovenous fistulas identification purposes and to confirm valve disruption and patency.

Storage Conditions

Store in a cool, dry place.

Restерilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

Safe Handling and Disposal

This device is a single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations:

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no heavy metals or radioisotopes, and is not infectious or pathogenic. Please consult local regulations to verify proper disposal of sharps.

Cleaning:

Devices considered necessary to return should be cleaned using one of the following:

- a. Sodium hypochlorite solution (500-600 mg/l), or
- b. Peracetic acid solution with subsequent ultrasonic treatment

Devices should then be decontaminated with either:

- a. 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
- b. Ethylene oxide gas

Devices should be completely dried prior to packaging.

Packaging:

Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of

maintaining the product without puncture of the packaging under normal conditions of transport.

The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.

Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.

Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.

Shipping paper and content marking for the outer shipping container are not required.

Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803, USA

Notices: Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd