# Syntel® Silicone Thrombectomy Catheter

Syntel® Silicone Thrombectomy Catheter
English — Instructions for Use
Cathéter de thrombectomie en silicone Syntel®

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### Syntel® Silicone Thrombectomy Catheter

(Model Numbers A4545, A4548, A4554, A4558, A4568)

English – Instructions for Use

# Rx only STERILE EO

Ref	Size	Maximum Liquid Capacity (ml)	CO2 Volume (ml)	Inflated Diameter (mm)	Recommended Syringe	
A4545	4F x 50cm	0.75	1.7	9	1cc	
A4548	4F x 80cm	0.75	1.7	9	1cc	
A4554	5F x 40cm	1.50	3.0	11	3сс	
A4558	5F x 80cm	1.50	3.0	11	3сс	
A4568	6F x 80cm	1.50	3.0	12	3сс	

### Description

The Syntel® Silicone Thrombectomy Catheter consists of a catheter body with a balloon on the distal end and an inflation hub on the proximal end.

### **Indications For Use**

The Syntel® Silicone Thrombectomy Catheter is indicated for use in vascular grafts and peripheral venous thrombectomy procedures.

### **Contraindications**

The Syntel® Silicone Thrombectomy Catheters are contraindicated for endarterectomy procedures or for use as a vessel dilator.

### **Precautions**

- Balloon degradation and rupture may result from deposits within the vessel, excessive handling or exposure to adverse environmental conditions
- Maximum recommended inflation volumes should not be exceeded

### **Complications:**

Complications associated with the use of vascular catheters include, but are not limited to: systemic infection, local hematomas, vessel dissection, perforation and rupture, bleeding, vessel thrombosis, distal emboli of blood clots or arteriosclerotic plaque, air embolus, aneurysms, vessel spasms, and balloon rupture or tip separation with fragmentation and distal embolization.

### Warnings

- Do not use if the product or sterile package is damaged
- Air should not be used to inflate the balloon

### **Instructions For Use**

- 1. The catheter balloon should be inflated with a sterile, nonparticulate, blood compatible liquid or CO<sub>2</sub>.
- 2. If using  $CO_2$ , corrections to inflation volumes may be necessary due to  $CO_2$  diffusion through the balloon wall.
- 3. Purge the balloon of all air with the inflation medium.
- 4. For optimal control, use the recommended syringe capable of holding the inflation volume. See Specification Table.
- 5. Insert the catheter through the arteriotomy/ venotomy, and position the balloon just beyond the occlusion.
- 6. Inflate the balloon with sterile medium. Inflation should continue until the balloon can be felt to have engaged the vessel wall.
- 7. Withdraw the catheter and thrombus through the arteriotomy/ venotomy.

### Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

### **How Supplied**

This device is supplied sterile and for single patient use. Pyrogen testing has been performed on this device. This device is supplied non-pyrogenic.

### **Storage Conditions**

Handle with care. Product should be stored in a clean, cool, dry area away from chemical fumes.

### **Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed of according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

### Cleaning

- 1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
- 2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
- 3. Devices should be completely dried prior to packaging.

### **Packaging**

- Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such
  packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product
  without puncture of the packaging under normal conditions of transport.
- 2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
- 3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
- 4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
- 5. Shipping paper and content marking for the outer shipping container are not required.
- 6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular Attn: Complaint Lab 63 Second Avenue Burlington, MA 01803, USA

### **Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular, Inc. Customer Service: Tel: (781) 221-2266
LeMaitre Vascular GmbH Otto-Volger-Str. 5a/b 65843 Sulzbach/Ts., Germany Tel: +49-(0)6196-659230
LeMaitre Vascular, Inc. 63 Second Avenue Burlington, MA 01803, USA