CardioCel® Bioscaffold Patch Instructions for Use

DESCRIPTION

CardioCel® is prepared from bovine pericardium.

The pericardium is procured from cattle originating in Australia. Australia is considered as being of negligible risk for BSE.

The device is a sterile, off-white, moist, pre-cut, flat or curved sheet of acellular collagen, presented sterile in a solution of propylene glycol and sealed in a container impermeable to air and moisture.

SIZES

The CardioCel bioscaffold patch comes in a range of sizes from 4cm² to 84cm² with curves that range from 30-90°. The catalogue numbers UCXXYY reflect the product dimensions where XX is the shortest dimension of the patch and YY is the longest dimension of the patch. The CardioCel Neo label is applied to products that are ≤ 0,4 mm thickness.

STERILIZATION

CardioCel is sterilized by exposure to, and storage in, a solution of 4% propylene oxide in sterile water (Osmolarity of ≤ 10 mOsm/L) which converts to propylene glycol solution.

INDICATIONS FOR USE

CardioCel is indicated for use in pericardial closure and the repair of cardiac and vascular defects including intracardiac defects; septal defects, valve and annulus repair; great vessel reconstruction, peripheral vascular reconstruction and suture line buttressing.

CONTRAINDICATIONS

CardioCel is not designed, sold or intended for use except as indicated.

The use of a tissue bioprosthesis must be considered on an individual patient basis.

INTENDED USER

The intended users of this device are qualified Cardiothoracic, Vascular and General Surgeons.

INTENDED PATIENT POPULATION

The CardioCel device is designed for permanent implantation in humans, and is indicated for repairing or reconstructing cardiac and vascular deformities.

WARNINGS

To avoid damage to the product, do not expose to any substances other than those specified in the rinse procedure.

Do not freeze.

Do not expose the device to extreme heat.

Do not re-sterilize.

Single use only and cannot be re-used.

Use of the device following a compromise in sterility may result in infection.

Do not use beyond the indicated expiration date.

Studies have not been performed to evaluate the safety or compatibility of this device during magnetic resonance imaging (MRI) scans. As such, the potential hazards of MRI testing on patients receiving this bio-prosthesis are unknown.

PRECAUTIONS

Store the package the right side up.

The outside of the jar is not sterile and must not be placed in the sterile field.

Do not use the device if the tamper evident seal is broken.

Do not use the device if the freeze indicator inside the lid of the device package indicates that the device has been exposed to temperatures of 0°C (32°F).

Do not use the device if there is evidence of damage to, or leakage from, the jar, as sterility of the product may have been compromised.

Do not use the device if the solution is turbid or coloured as this may be evidence that the product has been compromised.

Do not expose the patch to any solutions, chemicals, antibiotics, antimycotics, or other drugs except for the storage solution or sterile physiological saline, as irreparable damage to the patch may result that is not apparent under visual inspection.

Rinse surgical gloves to remove glove powder prior to touching CardioCel.

Prior to surgery, prospective patients or their representatives should be informed about possible complications which may be associated with the use of this device. Complete heart block has been reported for procedures involving cardiac repair near the atrial ventricular conduction bundles, most notably for repair of atrial septal defects.

As with any surgical procedures, infection is a possible complication. Monitor patient for infection and take appropriate therapeutic action.

PREPARATION FOR SURGICAL USE

Do not use if jar has evidence of having previously been opened or damaged.

First Handling

Remove jar from outer cardboard package and polystyrene inserts. Do not place the jar in the sterile field.

Inspect jar. Do not use if the tamper proof seal is missing or damaged, or there is evidence of moisture or leakage from the jar, or if the solution is turbid or coloured. The contents must be handled using aseptic techniques.

Opening

Open the jar without contaminating the contents. A member of the sterile team should remove the patch using sterile atraumatic forceps to grasp the edge of the patch and remove from the jar.

Immerse in sterile physiologic saline.

Rinsing

No extensive rinsing is required. Immersing the patch in sterile physiological saline is required to avoid dehydration of the patch prior to implantation.

THE DEVICE MUST REMAIN MOIST AT ALL TIMES.

Presentation to Surgeon and Use

The implant should be presented to the surgeon using sterile atraumatic forceps.

Surgical Use

The surgeon may cut and shape CardioCel to suit the requirements of the procedure.

Care should be taken when handling the device, for example by using atraumatic forceps, to avoid tearing or otherwise damaging the patch.

The device should be visually examined for damage noting it may have a smooth and a rough side.

CardioCel may be sutured or stapled in place. Being a strong material, it will take and hold sutures easily and firmly, and will remain in situ while it is incorporated into surrounding tissue.

When implanting by suture, suture bites should be taken 2 to 3 millimeters from the edge of the patch.

Note: After removal from its jar, all unused pieces of CardioCel should be discarded.

The solution in which the CardioCel is sterilized/stored can be disposed of according to hospital procedures for non-hazardous materials

STORAGE

CardioCel should be stored at room temperature, never below 2°C or above 25°C, and away from direct heat source. The implant should not be re-sterilized.

ADVERSE REACTIONS

The effect on long term surgical outcome of damaging bovine pericardium by chemicals or substances other than saline, by freezing, or by exposure to extreme heat, or chemical sterilization has not been investigated.

When bovine pericardium has been used for pericardial closure, cases of epicardial inflammatory reactions have been reported.

Other adverse events associated with bioprosthetic pericardial patches that have been reported in the literature have included: calcification; hemolysis; flow obstruction; thromboembolism; endocarditis; pericardial adhesions; inflammation; degeneration of the implants; and formation of clinically significant fibrous tissue.

LIMITED PRODUCT WARRANTY; LIMITATION OF REMEDIES

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES

NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device. IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction. A revision or issue date for these instructions is located in the bottom right corner of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

SYMBOL DEFINITIONS			
	Contents are Sterile		Product derived from Australian cattle.
	Single use only. Do not re-use.		Thickness
	Attention! See Instructions for Use.		Consult instructions for use.
	Minimum and maximum storage temperatures		Do not re-sterilize.
	Use by (with date)		Do not use package if product is damaged.
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.		Date of manufacture
	Catalogue Number		Manufacturer
	LOT Number		Serial Number.
	Sterilized using Propylene Oxide		

CardioCel and LeMaitre are registered trademarks of LeMaitre Vascular, Inc.

©2024 LeMaitre Vascular, Inc. All rights reserved.

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