**Pruitt® Aortic Occlusion Catheter**

English — Instructions for Use

**Pruitt® Aorta-Okklusjonskatheter**

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English – Instructions for Use

## Introduction

LeMaitre Vascular's Pruitt Aortic Occlusion Catheter offers the option to occlude the aorta via a transluminal approach when there is an acute aortic emergency, or when dissection and cross-clamping are difficult and time consuming (i.e. when the aneurysm is excessively large, dense, adherent to neighboring structures, such as the intestines, or especially if it is ruptured.) Intraluminal balloon occlusion may be accomplished by direct insertion through the wall of the aneurysm.

## Description

The Pruitt Aortic Occlusion Catheters are 12 French (4.0 mm), dual lumen catheters with a large, latex balloon (maximum liquid inflation capacity 50 mL) specifically designed and sized for use in the outlined general procedures. The first lumen (inflation lumen indicated by the white stopcock) is used for balloon inflation, while the second lumen (irrigation lumen indicated by the blue stopcock) allows access to the vessel distal to the occlusion. Other features include a stopcock with a luer-lock fitting at the proximal end of the irrigation lumen to facilitate control of such procedures and the inflation lumen to maintain balloon inflation level throughout the procedure, as well as a balloon wall thickness designed to reduce the possibility of puncture by calcium deposits.

A stainless steel stylet is inserted in the irrigation lumen of the catheter and serves as a stiffening medium to aid the physician during the introduction of the catheter into the patient's aorta.

The device is considered an orphan device in the European market and the premarket clinical data is relatively limited.

## Indication for Use

The Pruitt Aortic Occlusion Catheter is indicated to occlude the aorta to achieve control of blood flow during abdominal aortic vessel repair, aortic root replacement, and aorta arch repair procedures.

## Intended Use

The Pruitt Aortic Occlusion Catheter is intended to obtain rapid control of in-flow blood in the aorta in cases of ruptured aortic aneurysm or in other conditions when dissection of the neck of the aneurysm for different reasons may be especially difficult.

## Intended User

The Pruitt Aortic Occlusion Catheter is a surgical tool intended for use by experienced vascular surgeons trained in the procedures for which they are intended.

## Patient Population

Adults of any gender or ethnicity requiring treatment for aortic vessel repair, aortic root replacement, and aorta arch repair.

## Part of the Body Contacted

The Pruitt Aortic Occlusion Catheter will come in contact with the aorta.

## Clinical Condition

Aortic vessel repair, aortic root replacement, and aorta arch repair

## Clinical Benefits

The Pruitt Aortic Occlusion Catheter is used in the control of in-flow blood in the aorta, which indirectly benefits the patient by positively impacting their treatment management.

## Contraindications

1. The catheter is not to be used as a dilation catheter.
2. The catheter is not to be used for the introduction of drugs other than saline.
3. The catheter is a temporary device and cannot be implanted.

## Warnings

1. Do not reuse. The catheter is for single use only.
2. Air or gas should not be used to inflate the balloon during patient use.
3. Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED the recommended maximum balloon inflation capacity (maximum liquid inflation capacity 50 mL).
4. Exercise caution when encountering extremely diseased vessels. Arterial rupture or balloon failure due to sharp calcified plaque, may occur.
5. Deflate the balloon prior to inserting or withdrawing the catheter. Avoid using excessive force to push or pull catheter against resistance.
6. The possibility of balloon rupture or failure must be taken into account when considering the risk involved in a balloon catheterization procedure.
7. All agents to be infused should be used according to the manufacturer's Instructions for Use.
8. If the catheter is occluding blood flow to the kidneys, it should not be left in longer than 30-45 minutes.

## Precautions

1. Inspect the product and package prior to use and do not use the catheter if there is any evidence that the package or the catheter has been damaged.
2. Avoid extended or excessive exposure to fluorescent light, heat, sunlight, or chemical fumes to reduce balloon degradation. Excessive handling during insertion, or plaque and other deposits within the blood vessel may damage the balloon and can increase the possibility of balloon rupture.
3. Ensure proper connections between all syringes and hubs to avoid the introduction of air.
4. Do not grasp the balloon with instruments at any time to avoid damage to the latex.
5. Aspirate the irrigation lumen of the catheter during insertion until there is free back flow of blood from the catheter to reduce the chance of air embolism.

## List of Risks

*Potential device-related complications:*

- Aneurysms
- Aortic event
- Arterial Dissection
- Arteriovenous Fistula Formation
- Arterial Thrombosis
- Balloon Rupture

- Cerebrovascular accidents
- Distal Emboli of Blood Clots or Arteriosclerotic Plaque
- Hypertension or Hypotension
- Infection
- Intimal Disruption
- Mortality
- Paraplegia
- Post operative cardiac pulmonary arrest
- Renal Insufficiency
- Reoperation for dilated down stream
- Respiratory complications
- Tip Separation with Fragmentation and Distal Embolization
- Vessel Perforation and Rupture

*Potential procedure-related complications:*

- Air Embolus
- Arterial Spasm
- Hemorrhage
- Local Hematomas

**How Supplied**

The Pruitt Aortic Occlusion Catheter is supplied sterile and nonpyrogenic. This device is packaged in a non-sterile outer pouch and a sterile inner peel-open package. The sterility of the inner package is assured as long as it is unopened and undamaged.

**Procedure**

Pretest (perform before patient use)

1. Aspirate the balloon completely prior to inflation of the balloon with liquid.
2. Inflate the balloon with sterile saline and inspect for leaks. If there is any evidence of leaks around the balloon or if the balloon will not remain inflated, do not use the product.
3. Check the balloon by inflating and deflating with sterile saline for injection before use. If the balloon does not appear to function normally, do not use the product.

**General Use**

1. Flush the irrigation lumen (indicated by the blue stopcock) with sterile saline to remove air.
  2. With the balloon un-inflated introduce the catheter to the target anatomy over a 0.9mm "J" guidewire through an introducer sheath (7.3mm).
- NOTE:** The balloon catheter may be left in loco during graft placement and the proximal anastomosis carried out around the balloon catheter. If such a choice is made, thread the catheter through the graft BEFORE insertion.
3. During positioning, aspirate the irrigation lumen until there is free back flow of blood from the catheter, to reduce the chance of air embolism.
  4. Remove the guidewire and attach a 30 mL syringe loaded with sterile saline to the inflation lumen (indicated by the white stopcock).
  5. Open the stopcock and compress the syringe to inflate the balloon and occlude the vessel.

**WARNING:** Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED the recommended maximum balloon inflation capacity (maximum liquid inflation capacity 50 mL).

6. When occlusion is achieved, close the inflation stopcock to maintain balloon inflation.
7. To remove the device, open the white stopcock and use a 30ml syringe to draw saline from the balloon. Ensure device is fully deflated prior to removal.

**Storage Shelf Life**

The shelf life is indicated by the USE BY date on the package label. The Use By date printed on each label is NOT a sterility date. The Use By date is based on the normal life expectancy of the natural latex balloon when properly stored. The use of the catheter beyond the expiration date is not recommended because of potential balloon deterioration. LeMaitre Vascular, Inc. does not make provisions for replacing or reprocessing expired product.

Since natural rubber latex is affected by environmental conditions, proper storage procedures must be practiced to achieve optimum shelf life. The product should be stored in a cool dark area away from fluorescent lights, sunlight, and chemical fumes to prevent pre-mature deterioration of the rubber balloon. Proper stock rotation should be practiced.

**Re-sterilization/Repackaging**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

**Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

**Cleaning:**

1. Devices considered necessary to return should be cleaned using one of the following:
  - a. Sodium hypochlorite solution (500-600 mg/l), or
  - b. Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a. 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b. Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

**Packaging:**

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of

the primary receptacle. Cleaning methods should be detailed if possible.

3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803, USA

### **Summary of Safety and Clinical Performance**

To view the Pruitt Aortic Occlusion Catheter Summary of Safety and Clinical Performance document, please visit [www.lemaitre.com/sscp](http://www.lemaitre.com/sscp) then select the link "Pruitt Aortic Occlusion Catheter" to review the Pruitt Aortic Occlusion Catheter SSCP.

### **Notices: Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular GmbH
LeMaitre Vascular, Ltd.