XenoSure® Biologic Patch

English – Patient Information

XenoSure® Biologischer Patch

Deutsch – Patienteninformation

Patch biologique XenoSure®

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XenoSure® Biologic Patch

English – Patient Information

Biologic Vascular Patch Sizes	Width	Length	Model #
Biologic Vascular Patch - Tapered	0.6 cm	8 cm	0.6BV8M
Biologic Vascular Patch - Rounded	1 cm	6 cm	1BV6M
Biologic Vascular Patch - Tapered	0.8 cm	8 cm	0.8BV8M
Biologic Vascular Patch - Rounded	2 cm	9 cm	2BV9M
Biologic Vascular Patch - Tapered	1 cm	10 cm	1BV10M
Biologic Vascular Patch - Tapered	1.5 cm	10 cm	1.5BV10M
Biologic Vascular Patch - Tapered	1 cm	14 cm	1BV14M

Biologic Surgical Patch Sizes	Width	Length	Model #
Biologic Surgical Patch - Square	4 cm	4 cm	4BV4M
Biologic Surgical Patch - Rectangular	4 cm	6 cm	4BV6M
Biologic Surgical Patch - Rectangular	6 cm	8 cm	6BV8M
Biologic Surgical Patch - Rectangular	5 cm	10 cm	5BV10M
Biologic Surgical Patch - Rectangular	8 cm	14 cm	8BV14M
Biologic Surgical Patch - Rectangular	2.5 cm	15 cm	2.5BV15M
Biologic Surgical Patch - Rectangular	10 cm	16 cm	10BV16M

Intended Use

The XenoSure Biologic Patch is intended for use as a surgical patch material for heart repair, vessel reconstruction or vessel patching during surgical procedures.

Intended Patient Population

For narrowing neck arteries and weakened or damaged arteries in the thigh use: Adults of any gender or ethnicity who have narrowing neck arteries, and weakened or damaged arteries in the thighs.

For heart repair use: Pediatric patients who have holes in the upper or lower chambers of the heart requiring repair.

Self-care Instructions

1. Your new device is a foreign body and therefore needs close monitoring and careful observation. It may take 6-8 weeks for full healing.
2. After placement, the area may be swollen and tender for up to a week.
3. Watch for any new redness or tenderness.
4. Watch for any opening in the incisions.
5. Watch for numbness, tingling or pain near the new patch.

NOTE: If you experience any symptoms described in 3, 4 or 5 above please contact your provider.

6. Do not puncture or manipulate the patch.
7. If the patch was implanted in your leg, swelling in the extremity is expected because of increased blood flow. Elevate or move the extremity according to your provider's instructions.
8. Protect your patch from any trauma (sharp or blunt) objects.
9. It is preferable to have the surgical site covered for the first week to protect skin and incision(s). (Follow your provider's instructions.)
10. Keep bandages or wound covering on as per your provider's instructions.
11. If you have adhesive surgical tape or strips across your incision(s), wear loose clothing that does not rub against your incision(s). The adhesive surgical tape or strips will curl up and fall off on their own after a week.
12. You may shower or get the incision(s) wet, once your provider says you can. DO NOT soak, scrub, or have the shower beat directly on the incision(s).
13. DO NOT soak in the bathtub, a hot tub, or a swimming pool. Ask your provider when you can start doing these activities again.
14. Your provider will tell you how often to change your wound covering and when you may stop using one. Keep your incision(s) dry. If your incision(s) goes to your groin, keep a dry gauze pad over it to keep it dry.
15. Clean your incision(s) with soap and water every day once your provider says you can. Look carefully for any changes. Gently pat it dry.
16. DO NOT put any lotion, cream, or herbal remedy on your incision(s) without first discussing with your provider.
17. Consult your provider for instructions on taking any prescription or over-the-counter medications after your surgery.

Intended Performance

The patch is intended for heart or vessel reconstruction or repair.

Undesirable Side Effects

1. Electrical signal blockage within the heart can occur during repair procedures near cardiac nerve bundles.
2. Some chemicals (glutaraldehyde) may attack your immune system. The benefits of use of the XenoSure Biologic Patch have been weighed against the possible risk of late tissue weakening.
3. Animal studies with similar devices to the XenoSure patch have shown stiffening and signs of weakening over time. Findings include inflammation between your patch and surrounding tissue. This can lead to weakening of the patch.
4. Devices using similar material to your XenoSure Biologic Patch have been associated with inflammation and adhesions of the patch to the heart. The adhesions may increase the difficulty of repeat procedures.

Long-term Protection Measures for Your Patch

- For cardiac implants there are no long term protection measures. Please contact your provider for further guidance.
- If implanted in an extremity:
 - avoid prolonged extreme extension of the arm or leg with the implantation as it could lead to nerve damage.
 - avoid extreme or abrupt movements of the arm, shoulder, or legs during a post-operative period of 1.5 to 2 months. Specifically, you should not reach out in front, raise arms above shoulder level, throw, pull, stride, or twist.
 - avoid sleeping on the patch implantation side of your body as it may cause compression.
- If the patch is implanted on your neck, take care not to overextend your neck by bending or twisting frequently.

Risks Related to Interactions with Other Equipment

- None

Post-surgical Monitoring

- Check your incision(s) every day.
- Call your provider immediately if you have any signs of a blood clot, swelling, unusual skin color or infection, such as:
 1. Increased pain
 2. Swelling, redness, or red streaks
 3. Blood or pus draining from the access
 4. Numbness
 5. Fever
- If you have any symptoms of a stroke or neurological problems such as difficulty with speaking, your vision, or movement control, check with your provider immediately.

Lifetime of the Device

- The lifetime of the device has been shown to be safe and effective up to 12 years for the treatment of narrowing neck arteries, 5 years for the treatment of weakened or damaged arteries in the thighs, and 7 years in the treatment of heart repair. Proper care and regular medical follow up may help make this device last much longer.
- To ensure your patch functions as intended, follow the guidance of your health care provider.

When to Contact Your Provider

Closely watch for any changes in your health. Seek emergency care anytime you experience:

- Lost consciousness
- Trouble breathing
- Extremity has severe pain or becomes cold, pale, blue, tingly, or numb
- Pain that does not get better after you take pain medicine
- Loose stitches or your incision(s) opens
- Extensive bleeding from the incision(s)
- Signs of infection, such as:
 1. Increased pain, swelling, warmth, or redness
 2. Red streaks leading from the incision(s)
 3. Pus draining from the incision(s)
 4. A fever over 101°F (38.3°C)
 5. You are nauseous or cannot keep fluids down
- You have chest pain, dizziness, problems thinking clearly, or shortness of breath that does not go away when you rest
- You are coughing up blood or yellow or green mucus
- You have chills
- You experience abdominal pain or bloating

Device Materials

The patient may be exposed to the following materials and substances for worst-case largest patch 10 x 16 cm:

- Processed tissue from cow hearts up to 6.5 grams
- Low levels of residual glutaraldehyde up to 1.4 milligrams
- Low levels of residual formaldehyde up to 1.6 milligrams respectively
- Low levels of residual isopropyl alcohol up to 0.5 milligrams respectively

The XenoSure patch has passed testing to ensure it is safe.

Summary of Safety and Clinical Performance

To view the XenoSure Biologic Patch Summary of Safety and Clinical Performance document, please visit www.lemaitre.com/sscp then select the link "XenoSure Biologic Patch" to review the XenoSure SSCP.

In Case of an Emergency

Any serious adverse reaction should immediately be reported to your provider.

- Any device related emergency that occurs in the EU should be reported to the local Competent Authority of the country where you are located and the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or; <https://www.lemaitre.com/contact-us/email-us>).

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Dimensioni cerotto chirurgico biologico	Larghezza	Lunghezza	N. di modello
Biologic Surgical Patch - Square	4 cm	4 cm	4BV4M
Biologic Surgical Patch - Rectangular	4 cm	6 cm	4BV6M
Biologic Surgical Patch - Rectangular	6 cm	8 cm	6BV8M
Biologic Surgical Patch - Rectangular	5 cm	10 cm	5BV10M
Biologic Surgical Patch - Rectangular	8 cm	14 cm	8BV14M
Biologic Surgical Patch - Rectangular	2,5 cm	15 cm	2.5BV15M
Biologic Surgical Patch - Rectangular	10 cm	16 cm	10BV16M

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd