Artegraft™ Collagen Vascular Graft

Patient Information – English

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Intended Use

The Artegraft is intended to serve as a substitute channel for blood where replacement of blocked or diseased arteries is required or to establish a blood channel for getting dialysis.

Intended Patient Population

The product is designed for adult patients with variable weights, diagnoses, and health statuses.

Self-Care Instructions

1. Your new device is a foreign body and therefore needs close monitoring and careful observation. It may take 6-8 weeks for full healing.
2. After placement, the implant area may be swollen and tender for up to a week.
3. Watch for any new redness or tenderness.
4. Watch for any opening in the incisions.
5. Watch for numbness, tingling or pain in the leg, the new graft.

NOTE: If you experience any symptoms described in 3, 4, or 5 above, please contact your provider.

6. Do not puncture or manipulate the graft.
7. If the graft was implanted in your leg, swelling in the extremity is expected because of increased blood flow. Elevate or move the extremity according to your provider's instructions.
8. It is preferable to have the surgical site covered for the first week to protect skin and incision(s). (Follow your provider's instructions)
9. Keep bandages or wound covering on as per your provider's instructions.
10. If you have adhesive surgical tape or strips across your incision(s), wear loose clothing that does not rub against your incision(s). The adhesive surgical tape or strips will curl up and fall off on their own after a week.
11. You may shower or get the incision(s) wet, once your provider says you can. DO NOT soak, scrub, or have the shower beat directly on the incision(s).
12. DO NOT soak in the bath tub, a hot tub, or swimming pool. Ask your provider when you can start doing these activities again.
13. Your provider will tell you how often to change your wound covering and when you may stop using one. Keep your incision(s) dry. If your incision(s) goes to your groin, keep a dry gauze pad over it to keep it dry.
14. Clean your incision(s) with soap and water every day once your provider says you can. Look carefully for any changes. Gently pat it dry.
15. DO NOT put any lotion, cream, or herbal remedy on your wound without first discussing with your provider.
16. Consult your provider on taking any prescription or over-the-counter medications after your surgery.

Undesirable Side Effects

1. Electrical signal blockage within the heart can occur during repair procedures near cardiac nerve bundles.
2. Some chemicals (glutaraldehyde) may attack your immune system. The benefits of use of the Artegraft have been weighed against the possible risk of late tissue weakening.
3. Animal studies with similar devices to the Artegraft have shown stiffening and signs of weakening over time. Findings include inflammation between your patch and surrounding tissue. This can lead to weakening of the graft.
4. Devices using similar material to your Artegraft have been associated with inflammation and adhesions of the patch to the heart. The adhesions may increase the difficulty of repeat procedures.

Long-term Protection Measures for Your Graft

- For cardiac implants there are no long term protection measures. Please contact your provider for further guidance.
- Avoid prolonged extreme extension of the arm or leg with the implantation as it could lead to nerve damage.
- Avoid extreme or abrupt movements of the arm, shoulder, or legs during a post-operative period of 1.5 to 2 months. Specifically, you should not reach out in front, raise arms above shoulder level, throw, pull, stride, or twist.
- Avoid sleeping on the graft implantation side of your body as it may cause compression.

Risks Related to Interactions with Other Equipment

- None

Post-surgical Monitoring

- Check your incision(s) every day.
- Call your provider immediately if you have any signs of a blood clot, swelling, unusual skin color or infection such as:
 1. Increased pain
 2. Swelling, redness, or red streaks
 3. Blood or pus draining from the incisions
 4. Numbness
 5. Fever
- If you have any symptoms of a stroke or neurological problems such as difficulty with speaking, your vision, or movement control, check with your provider immediately.
- The nature and frequency of regular or preventive examination, monitoring or maintenance will be determined by your provider. This will be based on your underlying medical condition and the status of your graft.

Additional Post-surgical Monitoring for Hemodialysis

Check your access site for a pulse or “thrill” vibration. To feel it, place your first two fingers over the graft, or listen with a stethoscope. If you are able to listen with a stethoscope, you should be able to hear a repetitive whooshing sound. If you do not hear it or feel a pulse, the graft may be malfunctioning. Contact your provider with any concerns.

Lifetime of the Device

The Artegraft has demonstrated effectiveness for up to:

- 40 months for dialysis treatment
- 5 years for lower body blood vessel repair
- 32 months for the repair of damaged arteries

To ensure your graft functions as intended, follow the guidance of your health care provider.

When to Contact Your Provider

Closely watch for any changes in your health. Seek emergency care anytime you experience:

- Lost consciousness
- Trouble breathing
- Extremity has severe pain or becomes cold, pale, blue, tingly, or numb
- Pain that does not get better after you take pain medicine
- Loose stitches or your incision(s) opens
- Extensive bleeding from the incision(s)
- Signs of infection, such as:
 1. Increased pain, swelling, warmth, or redness
 2. Red streaks leading from the incision(s)
 3. Pus draining from the incision(s)
 4. A fever over 101°F (38.3°C)
 5. You are nauseous or cannot keep fluids down
- You have chest pain, dizziness, problems thinking clearly, or shortness of breath that does not go away when you rest
- You are coughing up blood or yellow or green mucus
- You have chills
- You experience abdominal pain or bloating

Device Materials

The patient may be exposed to the following materials and substances:

- Surgical sutures
- Blood vessel taken from cows
- Low levels of residual ethanol up to 0.36 g per graft
- Up to 3.6 mg of propylene oxide

The Artegraft Collagen Vascular Graft has passed testing to ensure its biocompatibility.

In Case of an Emergency

Any serious adverse reaction should immediately be reported to your provider.

- Any device related emergency that occurs in the EU should be reported to the Competent Authority of the country where you are located and the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or; <https://www.lemaitre.com/contact-us/email-us>).
- Any device related emergency that occurs in Australia should be reported to the Therapeutic Goods Administration via their website (<https://www.tga.gov.au/>) and the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or; <https://www.lemaitre.com/contact-us/email-us>).

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LeMaitre Vascular GmbH
LeMaitre Vascular, Ltd.