Artegraft™ Collagen Vascular Graft

English – Instructions for use

Artegraft™ Kollagen-Gefäßprothese

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Artegraft™ Collagen Vascular Graft

Artegraft™ Collagen Vascular Graft

(Model Numbers AG540M, AG630M, AG636M, AG640M, AG645M, AG715M, AG730M, AG735M, AG740M, AG745M, AG750M, AG830M, AG840M, AG845M, AG1015M, AG1030M)

English – Instructions for Use

Device Description

The Artegraft is composed of a section of specially selected bovine carotid artery to serve as a substitute conduit for blood where bypass or replacement of occluded or diseased arterial segments is required or to establish a conduit for hemodialysis.

Not Made with Natural Rubber Latex

Non-Antigenic

Indications For Use

The Artegraft is indicated for the following:

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Intended Use

The Artegraft is intended to serve as a substitute conduit for blood where bypass or replacement of occluded or diseased arterial segments is required or to establish a conduit for hemodialysis.

Intended User

The Artegraft is a surgical tool intended for use by experienced vascular surgeons trained in the procedures for which they are intended.

Patient Population

Adults of any gender or ethnicity that need an arterial bypass, an arteriovenous shunt, or surgery requiring an arterial graft.

Part of the Body Contacted

Arteries and veins in the lower and upper extremities.

Clinical Condition

Peripheral arterial disease, arteriovenous fistulae or any traumatic arterial injuries requiring bypass or interposition.

Clinical Benefits and Performance Characteristics

Clinical benefits of using the Artegraft for lower extremity bypass include enhanced graft anastomotic compliance due to vessel elasticity, reduced risk of graft infection, excellent long-term patency outcomes with a five-year overall primary patency rate of 66.7% and 74.6% secondary patency rate.

Clinical benefits of using the Artegraft for AV access include lower AV graft infection rates compared to ePTFE, access after 12 days, lower intervention rates as compared to ePTFE, and facilitates earlier removal of tunneled dialysis catheters.

The Artegraft is sourced from animal materials, as its performance is comparable to or better than grafts made from synthetic materials. The device's performance is also comparable to graft options made from other animal-derived (ovine or bovine) materials. It does not present any increased risk to the patient compared to alternatives, and its benefit-risk profile is favorable.

Contraindications

1. The Artegraft should not be used in venous or low pressure systems.
2. The Artegraft is also contraindicated in patients with known or suspected hypersensitivity to bovine collagen and bovine pericardium.

Lifetime of the Graft

The Artegraft has a lifetime of approximately:

- 40 months in the treatment of hemodialysis
- 5 years in the treatment of lower extremity bypass
- 32 months in the treatment of arterial trauma

Warnings

1. Do NOT re-sterilise the Artegraft. Use the graft immediately after opening the package and discard any unused portions.
2. The graft should not be used after the expiration date imprinted on the label.
3. After the Artegraft has been removed from the container in the manner prescribed to maintain its sterility, it should be gently and thoroughly washed and rinsed to minimize carryover of preserving fluid.
4. Silk is not recommended for anastomosis.
5. The graft is not to be used unless the capacity of the run-off vessel is adequate, as shown by pre-operative arteriography. Artegraft selection must be of comparable cross-sectional diameter to the host artery, particularly at the distal end, in order to avoid early thrombosis.
6. A minimum of twelve (12) days should be allowed after implantation before puncturing the graft with needles for hemodialysis. If edema appears around or distal to the graft, this should be allowed to resolve before cannulation.
7. Avoid external compression of the graft.

Precautions

1. In the event of early occlusion, re-exploration of the graft and removal of the thrombus with an embolectomy or thrombectomy catheter can be effective for restoration of long-term patency.
2. Patients subjected to heparin anticoagulant rinse and flush should be confirmed to be heparin-induced thrombocytopenia (HIT) free, and of heparin associated allergic reactions. Some surgeons recommend systemic heparinization of the patient after completion of the preparatory dissection, with or without subsequent neutralization with protamine sulfate. Others rely on the periodic injection of diluted heparin into the arterial tree during the period of vascular clamping and anastomosis. Post-operative heparinization is usually not employed.
3. During implantation ensure the graft does not twist as it is tunneled to its distal location.
4. Patients should be monitored for the presence of fever and transient low grade fever to prevent disruption of anastomosis
5. Patients should be monitored for high output congestive heart failure in patients with heart disease.

Potential Device-Related Complications:

- | | | |
|---------------------------|------------------|------------------------|
| • Aneurysm | • Infection | • Steal syndrome |
| • Bleeding | • Nonmaturation | • Stenosis |
| • Central venous stenosis | • Pseudoaneurysm | • Thrombosis formation |
| • Hematoma | • Seroma | |

Potential Procedure-related Complications

- Amputation
- Aneurysmal Degeneration
- Congestive Heart Failure
- Mortality
- Occlusion
- Reintervention

- Sepsis or Systemic Inflammatory Signs
- Wound Complication

How Supplied

The Artegraft is packaged in a sterile, single use, specially designed tube containing a sterilizing solution with 0.3% propylene oxide in 40% aqueous U.S.P. ethyl alcohol. Each tube is enclosed in a set-up box for protection during shipment and storage.

The sterility of the device is assured as long as the packaging is not opened or damaged. Do not use the device if the packaging has been opened or damaged. Once the package seal is broken, the graft should be used immediately. Any left-over material should be discarded.

Artegraft is not intended to be reused or resterilized. Attempts to reuse or resterilize the graft can result in contamination of the graft during the procedure.

The length and inner diameter of each Artegraft are specified on the packaging labels. The inner diameter of the Artegraft is approximate, rounded to the nearest millimeter (mm), due to the nature of the biologic source material. Product codes and sizes are referenced in the chart below.

NOTE: Outer diameters vary, but typically 1mm larger

Artegraft comes in the following sizes:

Model	Inner Diameter	Minimum Length	Model	Inner Diameter	Minimum Length
AG540M	4 mm	40 cm	AG740M	6 mm	40 cm
AG630M	5 mm	30 cm	AG745M	6 mm	45 cm
AG636M	5 mm	35 cm	AG750M	6 mm	50 cm
AG640M	5 mm	40 cm	AG830M	7 mm	30 cm
AG645M	5 mm	45 cm	AG840M	7 mm	40 cm
AG715M	6 mm	15 cm	AG845M	7 mm	45 cm
AG730M	6 mm	30 cm	AG1015M	8 mm	15 cm
AG735M	6 mm	35 cm	AG1030M	8 mm	30 cm

Storage

The graft must be kept in its original packaging until used. It must be kept in a dust free and dry environment at room temperature.

Removing the Graft from the container

1. To remove the Artegraft aseptically from its container, the container should be held in a semi-upright position.
2. After carefully removing the cap without contaminating the rim of the container, the supporting rod and the graft itself should be lifted from the container with long fine-tipped thumb forceps.
3. As soon as the graft itself becomes accessible, it should be grasped and removed from the container with the support rod and drain away excess alcohol from graft, care being taken not to touch the outside rim of the container.

Rinsing the Graft

1. Immerse the graft in a large basin with a minimum of 500 ml of sterile 0.9% sodium chloride solution (sterile saline).
2. Irrigate the graft's lumen at least seven times with ~30ml of solution using a bulb or piston syringe (with slip tip).
3. A final irrigation and immersion should be performed. In a second basin, prepare a suggested ratio of 120ml sterile saline solution with 2500 units of heparin, unless HIT is diagnosed (see PRECAUTIONS).
4. Before implanting, occlude one end of the graft and pressure test with a bulb or syringe filled with saline.

The graft is now considered ready for insertion into a patient.

Implantation

1. Accepted vascular surgical techniques should always be employed. Precise anastomosis without tension, bulging or twisting is mandatory.
2. The use of fine non-cutting rather than cutting-type needles helps avoid troublesome punctate bleeding at the anastomotic sites, particularly in the smaller diameter grafts.
3. If troublesome bleeding does occur, it may readily be controlled by simply applying a small piece of excess Artegraft as a patch to the bleeding area in a reverse manner so that the adventitia of the patch is in direct contact with adventitia.
4. In order to avoid possible snarling of the suture in the rather shaggy adventitia of the graft, it has been found convenient to start the stitch in the host artery rather than starting it in the graft. In so doing, care should be exercised to avoid the possibility of displacing an atheromatous plaque.
5. The graft lot number and its length and diameter should be recorded on the patient's chart, as well as the actual length of the graft implanted into the patient.

Patient Implant Resources

The Artegraft Collagen Vascular Graft is supplied with a Patient Implant Card (PIC). Please supply the patient with the completed PIC (instructions below) after implantation:

1. The front of the PIC supplied is to be completed by the operating surgeon/team.
2. There are 3 lines of information to be completed. Line #1 is for patient identification (e.g. patient name). Line #2 is for the operation date. Line #3 is for the address of the health care center or doctor where medical information about the patient may be found.
3. The back of the PIC contains the product and Manufacturer information.

After completion, provide the PIC to the patient for future reference.

Device Materials

The patient may be exposed to the following:

1. Device materials for worst-case largest graft 7 mm x 45 cm:
 - Acellular bovine carotid artery and Non-Absorbable Surgical Suture: 12.38 g per graft
2. Chemical residues for worst-case largest graft 7 mm x 45 cm:
 - Ethanol: up to 0.36 g per graft
 - Propylene oxide: up to 3.6 mg per graft

Safe Handling and Disposal

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located. This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. There are no special requirements for disposal. Please consult local regulations to verify proper disposal.

Dispose of storage solution according to local and federal regulations. Solution should not be disposed using septic systems. If there are no disposal restrictions, solution may be diluted and disposed in a sanitary sewer system.

Packaging and shipping of explanted Artegraft:

Return of the shipment to LeMaitre Vascular depends on 3 crucial questions:

1. Is the explant obtained from a patient with a known or presumed pathogenic condition at the time of the explant?
2. Is the explant obtained from a patient that has a known treatment history that involves therapeutic radionuclides within the last 6 months?
3. Has the clinician obtained consent from the patient for the sample to be returned to the manufacturer for research purposes?

If the response to question 1 or 2 is affirmative, LeMaitre Vascular does not provide adequate guidance for shipment. THESE EXPLANTS SHOULD NOT BE RETURNED TO LEMAIRE VASCULAR UNDER ANY CIRCUMSTANCES. In these cases, the explant should be disposed according to local regulations.

For explants that do not have pathogenic or radiologic hazards, please use the following:

Pre-explantation:

1. If possible, perform a CT or Ultrasound scan of the device to document patency.
2. LeMaitre Vascular can accept clinical information that is patient-anonymized. LeMaitre Vascular requests information including:
 - a) The original diagnosis which resulted in use of the implant.
 - b) The patient's medical history relevant to the implant, including the hospital or clinic at which the device was implanted.
 - c) The patient's implant experience prior to implant removal.
 - d) The hospital or clinic at which the explantation was performed, and date of retrieval.

Explanation:

1. Explanted Artegrafts should be directly transferred to a sealable container filled with a solution of **alkaline buffered 2% glutaraldehyde or 4% formaldehyde** prior to shipment.
2. Cleaning of explanted grafts should be minimal if necessary. Proteolytic digestion should not be used in any circumstances.
3. Artegraft explants should not be decontaminated under any circumstances. DO NOT autoclave the sample or use ethylene oxide gas to decontaminate.

Packaging:

1. Explants should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. Material that is absorbent and is cushioning should be selected for isolating the sealable container inside the secondary packaging. Primary and secondary packaging must then be packaged inside an outer package.
2. Explants in sealed primary containers should be labeled with an ISO 7000-0659 Biohazard symbol. The same symbol should be attached to the secondary packaging and to the outer packaging. Outer packaging should also be labelled with Name, Address and Telephone Number of Sender, and the statement, "Upon discovery of damage or leakage, the package should be isolated and sender notified".
3. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803, USA

Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. The unused section of the Artegraft prosthesis or the explanted prosthesis should be considered as biohazard and be disposed per hospital policy or local regulations.

Summary of Safety and Clinical Performance

To view the Artegraft Summary of Safety and Clinical Performance document, please visit www.lemaitre.com/sscp then select the link "Artegraft Collagen Vascular Graft" to review the Artegraft SSCP.

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAIRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAIRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAIRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAIRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAIRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

Artegraft™ Kollagen-Gefäßprothese

(Modellnummern AG540M, AG630M, AG636M, AG640M, AG645M, AG715M, AG730M, AG735M, AG740M, AG745M, AG750M, AG830M, AG840M, AG845M, AG1015M, AG1030M)

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd