AlboGraft® Polyester Vascular Graft

English – Patient Information

AlboGraft® Polyester-Gefäßprothese

Deutsch – Patienteninformationen

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English – Patient Information

Intended Use

AlboGraft Vascular Grafts are intended to be used as a replacement for diseased blood vessels.

Intended Patient Population

Adults of any gender or ethnicity in need of repair or replacement of diseased blood vessels or during surgery to correct blood flow.

Care Instructions

1. Your new device is a foreign body and therefore needs close monitoring and careful observation. It may take 6-8 weeks for full recovery.
2. After placement, the implant area maybe swollen and tender for up to a week.
3. Watch for any new redness or tenderness.
4. Watch for any opening in the incisions.
5. Watch for numbness, tingling or pain in the leg.

NOTE: If you experience any symptoms described in 3,4, or 5 above please contact your provider.

6. Do not puncture or manipulate the site of the incision or graft, including if visible, without medical staff supervision.
7. Swelling in the extremity is expected because of increased blood flow. Move according to your provider's instructions. Temporarily keep your leg elevated above your heart until swelling has gone down. Contact your care provider with any questions or concerns.
8. It is preferable to have the graft covered for the first week to protect skin and incisions. (Follow your provider's instructions).
9. Keep bandages or compression bandages on as per your provider's instructions.
10. If your staples have been removed, you will probably have surgical tape or strips across your incision. Wear loose clothing that does not rub against your incision. The adhesive surgical tape or strips will curl up and fall off on their own after a week.
11. You may shower or get the incision wet, once your provider says you can. DO NOT soak, scrub, or have the shower beat directly on them.
12. DO NOT soak in the bathtub, a hot tub, or swimming pool. Ask your provider when you can start doing these activities again.
13. Your provider will tell you how often to change your wound covering and when you may stop using one. Keep your incisions dry. If your incision goes to your groin, keep a dry gauze pad over it to keep it dry.
14. Clean your incision with soap and water every day once your provider says you can. Look carefully for any changes. Gently pat it dry.
15. DO NOT put any lotion, cream, or herbal remedy on your wound without first discussing it with your provider.

Long-term protection measures for your graft

- Avoid protracted extreme extension of the arm. Prolonged extreme extension may lead to nerve damage.
- Avoid extreme or abrupt movements of the arm, shoulder, or legs during a post-operative period of 1.5 to 2 months. Specifically, you should not reach out in front, raise arms above shoulder level, throw, pull, stride, or twist.
- Avoid sleeping on the graft implantation side of your body or crossing your legs for prolonged periods as it may cause compression.

Risks Related to Interactions with Other Equipment

None.

Post-surgical monitoring

- Check your incision every day.
- Call your provider immediately if you have any signs of a blood clot, swelling, unusual skin color, or infection, such as:
 1. Increased pain
 2. Swelling, redness, or red streaks
 3. Blood or pus draining from the incisions
 4. Numbness
 5. Fever
- Check for signs of good circulation. Your foot or leg should not be cool, pale, experience pain, or have other symptoms similar to pre-surgery. Call your provider if you are experiencing any of these symptoms.
- The nature and frequency of regular or preventive examination, monitoring or maintenance will determined by your provider. This will be based on your underlying medical condition and the status of your graft.

Lifetime of the Device

- The lifetime of the device has been shown to be safe and effective for up to 10 years. Proper care and regular medical follow-up may extend the life of the device.
- To ensure your graft functions as intended follow the guidance of your health care provider.

When to contact your Provider

Closely watch for any changes in your health. Seek emergency care anytime you experience:

- Lost Consciousness
- Trouble breathing
- Extremity has severe pain or becomes cold, pale, blue, tingly, or numb
- Pain that does not get better after you take pain medicine
- Loose stitches, or your incision opens
- Extensive bleeding from the incisions
- Signs of infection, such as:
 - a. Increased pain, swelling, warmth, or redness
 - b. Red streaks leading from the incision
 - c. Pus draining from the incision

- d. A fever over 101°F (38.3°C)
- e. You are nauseous or cannot keep fluids down
- You have chest pain, dizziness, problems thinking clearly, or shortness of breath that does not go away when you rest
- You are coughing up blood or yellow or green mucus
- You have chills
- You experience abdominal pain or bloating

Device Materials

The following materials are what make up the AlboGraft Polyester Vascular Graft:

1. Cowhide up to 9.0 grams
2. Glycerin up to 35.0 grams
3. Polyester Yarn up to 36.0 grams
4. Polypropylene up to 3.0 grams (ASC grafts only)
5. Polyethylene up to 2.0 grams (ASC grafts only)

The AlboGraft Polyester Vascular Graft has passed testing to ensure it is safe.

In case of an emergency

- Any serious adverse reaction should immediately be reported to your provider.
- Please make sure any serious incident that occurs in relation to the graft is reported to the Competent Authority of the country where you are located
- You may also report the incident to the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or <https://www.lemaitre.com/contact-us/email-us>)

Summary of Safety and Clinical Performance

To view the AlboGraft® Polyester Vascular Graft Summary of Safety and Clinical Performance document, please visit: www.lemaitre.com/sscp then select the link "AlboGraft® Polyester Vascular Graft" to review the AlboGraft SSCP.

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LeMaitre Vascular GmbH
LeMaitre Vascular, Ltd.