VascuTape® Radiopaque Tape

English - Instructions for Use

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English – Instructions for Use

Description:

VascuTape® Radiopaque Tape is a sterile medical device that can be used during vascular procedures (eg, stent implantation, percutaneous transluminal angioplasty, endovascular atherectomy, vena cava filter placement, stent-graft bypass, in situ bypass, etc.) or in other procedures that require x-ray positioning. VascuTape® Radiopaque Tape is for use on intact skin only. The markings on the radiopaque tape are for reference only, and are not calibrated.

Indication for Use:

The VascuTape Radiopaque Tape is indicated as a reference for measurements during imaging, including angiography.

Intended Use:

VascuTape Radiopaque Tape is intended to be placed on the patient's skin to assist during imaging procedures by serving as a frame of reference only. This product is used for reference measurement of objects in an x-ray image, and/or to understand the location of objects in an image in relation to where the tape is adhered to the patient's body.

Intended User:

The intended users of VascuTape Radiopaque Tape include vascular, general, and orthopedic surgeons, interventional radiologists, and interventional cardiologists.

Patient Population:

Adults of any gender or ethnicity requiring imaging such as angiography.

Part of the Body Contacted:

Patient skin; topical.

Clinical condition:

Patients needing imaging.

Clinical Benefits:

The VascuTape Radiopaque Tape is anticipated to have a beneficial effect on the health of the patient, resulting in meaningful (relevant to the patient) and measurable positive impacts on clinical outcomes by serving as a frame of reference during radiological procedures. This benefit is indirect.

Contraindications:

- Radiopaque tape should not be used with patients who are allergic to any of the materials that are in the non-tackified Acrylate adhesive.
- Radiopaque tape should not be implanted.

Warnings:

1. Do not use if the package is opened or damaged.
2. The product should be removed within 12 hours of application to the patient.
3. Do not reuse. Radiopaque Tape is for single use only.

Potential Complications:

As with all devices and surgical procedures, complications may occur. These include but are not limited to skin irritation and injury to the patient.

Procedure:

1. The radiopaque tape is packaged sterile. Take the radiopaque tape out of the packaging and discard the back sticker to expose the adhesive. (Fig. 1)
2. Apply the radiopaque tape to the patient's skin. (Fig. 2)
3. The radiopaque tape serves as an internal reference during the X-ray examination. (Fig. 3)
4. After the procedure is completed, or within a maximum of 12 hours of application to the patient, the product should be removed. (Fig. 4)

Re-sterilization/Repackaging:

This device is single-use only. Do not reuse, reprocess, or re-sterilize.


The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

Incident Reporting:

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

Safe Handling and Disposal:

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
 - a. Sodium hypochlorite solution (500-600 mg/l), or
 - b. Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
 - a. 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
 - b. Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803, USA

Limited Product Warranty; Limitation of Remedies:

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular, Ltd.
LeMaitre Vascular GmbH