Pruitt® Occlusion Catheter

English – Instructions for Use

Pruitt® Okklusionskatheter

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Cathéter d'occlusion Pruitt®

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Catetere per occlusione Pruitt®

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Catéter de oclusión Pruitt®

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English – Instructions for Use

Description

The Pruitt® Occlusion Catheter (POC) is specifically designed for the temporary occlusion of a blood vessel without the use of a clamp or ties. The latex balloon conforms to the diseased interior wall of the vessel providing occlusion with less damage to the vessel wall and lining. The Pruitt Occlusion Catheter is offered in 3F, 4F, and 5F sizes. A stopcock is bonded to the proximal arms of the Pellethane catheter to allow the balloon to remain inflated if desired. Depth markings of 50mm, 100mm, and 150mm on the distal end of the catheter aid the physician in positioning the balloon in the vessel.

Storage

The Occlusion Catheters should be stored in a cool dark area away from fluorescent lights and sunlight to prevent premature deterioration of the latex balloon. The recommended maximum shelf life (use by date) is shown on the label.

Indication for Use

The Pruitt Occlusion Catheter is indicated to occlude peripheral arteries to facilitate the performance of the thoracoabdominal aortic aneurysm (TAA) procedure.

Intended Use

The POC is intended to conform to the diseased interior wall of smaller vessels, including tibial, popliteal, lumbar and celiac blood vessels providing occlusion during open vascular repair procedures to control bleeding.

Intended User

The POC is a surgical tool intended for use by experienced vascular surgeons trained in the procedures for which they are intended.

Patient Population

Patients of any gender, age or ethnicity undergoing repair of aneurysmal and occlusive diseases requiring open vessel repair.

Part of the body contacted

Peripheral vessels (e.g. renal, mesenteric and other vessels of appropriate size)

Clinical Condition

Aneurysmal and occlusive disease

Clinical Benefits

The catheters are designed to assist the surgeon by providing temporary occlusion while the vascular repair procedures are being performed. As such, the clinical benefit provided by the devices is assessed primarily by evaluating whether the device performed as intended (i.e., successfully occluded blood flow intraoperatively).

Contraindications

1. The catheter is not to be used as a dilation catheter.
2. The catheter is a temporary device and cannot be implanted.

Warnings

1. The device is for single use only. Do not reuse, reprocess, or resterilize.
2. Air or gas should not be used to inflate the balloon if there is a possibility of embolization with balloon rupture.
3. Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED the recommended maximum balloon inflation capacity (See Specifications Table 1).
4. Deflate the balloon prior to withdrawing the catheter.
5. The possibility of balloon rupture or failure must be taken into account when considering the risk involved in a balloon catheterization procedure.
6. The maximum inflation cycle of this device is 10 cycles.
7. Do not exceed the maximum inflation time of 6 hours.

Precautions

1. Inspect the product and package prior to use and do not use the catheter if there is any evidence that the package or the catheter has been damaged.
2. Avoid extended or excessive exposure to fluorescent light, heat, sunlight, or chemical fumes to reduce balloon degradation.
3. Ensure proper connections between all syringes and hubs to avoid the introduction of air.
4. To avoid damage, do not grasp the balloon with instruments during insertion.
5. Aspirate before inflation.
6. Pretest the catheter before use: a) inflate the balloon to the recommended capacity with air and immerse the balloon in sterile water. If there is any evidence of air escaping around the balloon or if the balloon will not remain inflated, do not use the product. b) Also, check the balloon integrity by inflating and deflating with sterile saline for injection before use. If the balloon does not appear to function normally, do not use the product.
7. Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Instructions for Use

1. The POC is supplied in sterile peel-open packages. Inspect the package and do not use the catheter if there is any evidence that the package has been punctured or that the catheter has been damaged.
2. Inflate the balloon to the recommended capacity with air and immerse the balloon in sterile water. If there is any evidence of air bubbles escaping around the balloon or if the balloon will not remain inflated, do not use the catheter.
3. Aspirate the catheter with proper volume of sterile saline.
4. Insert the catheter into the target vessel while the balloon is deflated.
5. The depth markings indicate how far the catheter has been introduced into the vessel. If the catheter is at the 50mm depth marking, this indicates the catheter has been inserted 50mm into the vessel.
6. At the desired location, inflate the balloon until full occlusion is achieved. Do not over-inflate the balloon beyond the maximum volume listed in the Specification Table.
7. Close the stopcock to maintain the inflation of the balloon during the procedure.
8. After the procedure is completed, deflate the balloon by opening the stopcock and allow the saline to drain.
9. Remove the catheter from the patient and dispose according to the hospital's policy.

Potential Complications

- Air embolus
- Aneurysm
- Arterial dissection
- Arterial spasm
- Arterial thrombosis
- Arteriovenous fistula formation
- Balloon rupture
- Death
- Embolism
- Embolization of blood clots, arteriosclerotic plaque, or air
- Hemorrhage
- Hypertension or hypotension
- Infection
- Intimal disruption
- Local hematoma
- Paraplegia
- Renal insufficiency
- Stroke

- Tip separation with fragmentation and distal embolization

- Vessel perforation and rupture

Specifications			
Model	2103-36M	2103-46M	2103-56M
Maximum Liquid Capacity	0.5 mL	0.5 mL	1.0 mL
Inflated Diameter	7mm	9mm	11mm
Usable Length	27 cm		
Contents	One Catheter, One Syringe, 3 mL		

Re-sterilization/Repackaging

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

Incident Reporting

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
 - a. Sodium hypochlorite solution (500-600 mg/l), or
 - b. Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
 - a. 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
 - b. Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd.