## **Distal Perfusion Catheter**

English — Instructions for Use

## **Katheter zur distalen Perfusion**

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## **Cathéter de Perfusion Distale**

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**STERILE EO Rx only** LATEX

## **Product Description**

The Distal Perfusion Catheter is a 12 French (4.0 mm) catheter with a large through lumen and a distal latex balloon. The balloon inflation lumen, indicated by the white stopcock, is used for inflating the balloon, while the second lumen, indicated by the removable blue stopcock, allows for the flow of blood to the recipient vessel. The stopcock on the through lumen is removable to provide an option for the user to use the stopcock to provide flow control or remove the stopcock to improve the flow rate.

## **How Supplied**

The Distal Perfusion Catheter is supplied sterile and nonpyrogenic. This device is packaged in a non-sterile outer pouch and a sterile inner peel-open package. The sterility of the inner package is assured as long as it is unopened and undamaged.

## **Indication**

The Distal Perfusion Catheter is indicated for use as a distal perfusion conduit during open abdominal and thoracoabdominal aortic aneurysm repair.

## **Contraindications**

1. The catheter is not to be used for embolectomy, thrombectomy or vessel dilatation.
2. The catheter is a temporary device and should not be implanted.
3. The catheter is not to be used for the introduction of drugs other than saline, heparin, or contrast media.

## **Warnings**

1. The device is for single use only. Do not reuse, reprocess, or resterilize.
2. Air or gas should not be used to inflate the balloon during patient use.
3. Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED THE RECOMMENDED MAXIMUM BALLOON INFLATION CAPACITY (See Table 1).
4. Exercise caution when encountering extremely diseased vessels. Arterial rupture or balloon failure due to sharp calcified plaque may occur.
5. Deflate the balloon prior to inserting or withdrawing the catheter. Avoid using excessive force to push or pull catheter against resistance.
6. The possibility of balloon rupture or failure must be taken into account when considering the risk involved in a balloon catheterization procedure.
7. All agents to be infused should be used according to the manufacturer's Instructions for Use.

## **Precautions**

1. Inspect the product and package prior to use and do not use the catheter if there is any evidence that the package or the catheter has been damaged.
2. Avoid extended or excessive exposure to fluorescent light, heat, sunlight, or chemical fumes to reduce balloon degradation.
3. Aspirate the balloon prior to inflation to minimize the amount of air injected into the balloon.
4. Excessive handling during insertion, or plaque and other deposits within the blood vessel, may damage the balloon and can increase the possibility of balloon rupture.
5. Ensure proper connections between all syringes and hubs to avoid the introduction of air.
6. Do not grasp the balloon with instruments at any time to avoid damage to the latex.
7. Flush the perfusion lumen of the catheter with either saline or blood prior to insertion of the distal end to reduce the chance of air embolism.

## **Potential Complications**

As with all catheterization and surgical procedures, complications may occur. These may include, but are not limited to:

- |   |   |
|---|---|
| • Embolization of blood clots, arteriosclerotic plaque or air | • Aneurysms   |
| • Hypertension or hypotension                                 | • Arterial spasms   |
| • Infection   | • Local hematomas   |
| • Intimal disruption  | • Arteriovenous fistula formation                           |
| • Arterial dissection   | • Balloon rupture   |
| • Vessel perforation and rupture                              | • Tip separation with fragmentation and distal embolization |
| • Hemorrhage  | • Paraplegia  |
| • Arterial thrombosis   |   |

## **Procedure**

### **Pretest (perform before patient use)**

1. Aspirate the balloon completely prior to inflation of the balloon with liquid.
2. Inflate the balloon with sterile saline and inspect for leaks. If there is any evidence of leaks around the balloon or if the balloon will not remain inflated, do not use the product.
3. Check the balloon by inflating and deflating with sterile saline for injection before use. If the balloon does not appear to function normally, do not use the product.

## **General Use**

*Note: A removable blue stopcock is provided with the Distal Perfusion Catheter. This stopcock is attached to the female luer at the proximal end of the blood flow lumen. The stopcock is provided with the handle of the rotating center portion lined up with the flow path stopcock (the open position). To stop the flow, turn the handle 90 degrees in either direction so that the handle is perpendicular to the flow path of the stopcock.*

1. Attach the proximal (stopcock/luer end) of the perfusion lumen to the distal end of the blood flow circuit. If you have not primed the perfusion lumen with saline, be sure to flush the lumen with blood to eliminate the risk of air embolus when the catheter is secured in the vessel and blood flow is initiated.
2. Fill the 5 ml syringe with the appropriate amount of sterile saline (see Table 1) and attach it to the white stopcock of the balloon inflation lumen.
3. Insert the distal tip of the catheter with the balloon into the recipient vessel and inflate the balloon using the least volume necessary to secure the catheter in the vessel, taking care not to over inflate the balloon. Do not exceed the recommended maximum liquid capacity (see Table 1). Verify that the balloon is secure in the vessel before starting blood flow.

**Table 1: Recommended Maximum Liquid Capacity**

Model Number	Maximum Liquid Capacity	Inflated Diameter
2105-15	3.4 ml	20mm

4. After establishing blood flow through the Distal Perfusion Catheter, verify there is no leakage around the balloon.

5. Prior to removal of the catheter, first verify the blood flow has been stopped. Then, deflate the balloon by opening the white stopcock and drawing a vacuum with the syringe. Be sure the balloon is completely deflated before withdrawing the distal tip of the catheter from the vessel.

#### **Table 2: Flow Rates**

LeMaitre Vascular performed flow rate tests to simulate blood flow through the catheter at 120/80mmHg.

	Flow Rate without Removable Stopcock	Flow Rate with Removable Stopcock
Maximum	294 ml/min	279 ml/min
Minimum	282 ml/min	265 ml/min

#### **Storage Shelf Life**

The shelf life is indicated by the Use By date on the package label. The Use By date printed on each label is not a sterility date. The Use By date is based on the normal life expectancy of the natural latex balloon when properly stored. The use of the catheter beyond the expiration date is not recommended because of potential balloon deterioration. LeMaitre Vascular, Inc. does not make provisions for replacing or reprocessing expired product. Since natural rubber latex is affected by environmental conditions, proper storage procedures must be practiced to achieve optimum shelf life. The product should be stored in a cool dark area away from fluorescent lights, sunlight, and chemical fumes to prevent pre-mature deterioration of the rubber balloon. Proper stock rotation should be practiced.

#### **Resterilization/Re-use**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

#### **Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed accorded to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

#### **Cleaning:**

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

#### **Packaging:**

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803, USA

#### **References**

1. Estrera AL, et al., Neurologic outcome after thoracic and thoracoabdominal aortic aneurysm repair. Annals of Thoracic Surgery, 72:1225-1231, 2001.

#### **Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre is a registered trademark of LeMaitre Vascular, Inc.
©2025 LeMaitre Vascular, Inc. All rights reserved.
LeMaitre Vascular GmbH
LeMaitre Vascular Ltd