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Standardin nimi	Standardin viite: tarkistusvuosi
Sterilization of medical devices. Requirements for medical devices to be designated "STERILE". Part 2: Requirements for aseptically processed medical devices	EN 556-2:2015
Information supplied by the manufacturer of medical devices	EN 1041:2008
Cardiovascular implants and extracorporeal systems – Vascular prostheses -- Tubular vascular grafts and vascular patches	ISO 7198:2016
Biological evaluation of medical devices – Part 1: Evaluation and testing	ISO 10993-1:2009
Biological evaluation of medical devices – Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	ISO 10993-3:2009
Biological evaluation of medical devices – Part 4: Selection of tests for interactions with blood	EN ISO 10993-4:2006
Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	ISO 10993-5:2009
Biological evaluation of medical devices – Part 6: Tests for local effects after implantation	EN ISO 10993-6:2007
Biological evaluation of medical devices – Part 10: Tests for irritation and delayed-type hypersensitivity	ISO 10993-10:2010
Biological evaluation of medical devices – Part 11: Tests for systemic toxicity	ISO 10993-11:2018
Biological evaluation of medical devices Part 17: Establishment of allowable limits for leachable substances	EN ISO 10993-17:2008
Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	ISO 11607-1:2006
Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes	ISO 11607-2:2006
Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products	ISO 11737-1:2006
Tests of sterility performed in the definition, validation and maintenance of a sterilization process	ISO 11737-2:2009
Aseptic processing of health care products – Part 1: General requirements	ISO 13408-1:2008
Lääkinnälliset laitteet – Laadunhallintajärjestelmät – Vaatimukset viranomaismääräyksiä varten	EN ISO 13485:2016

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Viitteet:

1. Chuatrakoon B, Nantakool S, Rerkasem A, Orrapin S, Howard DPJ, Rerkasem K. Routine or selective carotid artery shunting for carotid endarterectomy (and different methods of monitoring in selective shunting). *Cochrane Database of Systematic Reviews*. 2022;2022(6).
2. Vuurberg NE, Post IC, Keller BPJ, Schaafsma A, Vos CG. A systematic review and meta-analysis on perioperative cerebral and hemodynamic monitoring methods during carotid endarterectomy. *Annals of Vascular Surgery*. 2023;88:385-409.
3. Cho JS, Song S, Huh U, et al. Comparing carotid endarterectomy and carotid artery stenting: retrospective single-center analysis. *Ann Palliat Med*. 2022;11:3409-3416.
4. Kazantsev A, Korotkikh A, Lider R, et al. Results of carotid endarterectomy with the use of temporary shunts with reduced retrograde pressure in the internal carotid artery—analysis of the multicenter Russian register. *Indian Journal of Thoracic and Cardiovascular Surgery*. 2023;39(3):244-250.
5. Levin SR, King EG, Farber A, Cheng TW, Rybin D, Siracuse JJ. Unplanned Shunting Is Associated with Higher Stroke Risk after Eversion Carotid Endarterectomy. *Annals of Vascular Surgery*. 2022;87:362-368.
6. Ribieras AJ, Tabbara M, Rey J, Velazquez OC, Bornak A. Outcomes and role of shunting during carotid endarterectomy for symptomatic patients. *Journal of vascular surgery*. 2022.
7. Soliman SS, Lazar A, Millner NB, et al. Comparing carotid endarterectomies with or without shunting in symptomatic and asymptomatic patients. *The American Journal of Surgery*. 2023.
8. Squizzato F, Siracuse JJ, Shuja F, et al. Impact of Shunting Practice Patterns during Carotid Endarterectomy for Symptomatic Carotid Stenosis. *Stroke*. 2022;53.
9. Ucci A, de Troia A, D'Ospina RM, et al. Carotid endarterectomy in asymptomatic octogenarians: Outcomes at 30 days and 5 years. *Vascular*. 2023;31(1):98-106.
10. Yuan W, Huo R, Ma K, et al. A single-center retrospective study with 1-year follow-up after CEA in patients with severe carotid stenosis with contralateral carotid artery occlusion. *Frontiers in Neurology*. 2022;13:971673.
11. Zhao W, Gao F, Wu C, et al. Severe contralateral carotid stenosis or occlusion drive 30-day risk after carotid endarterectomy. *Vascular*. 2022;30(1):3-13.
12. Aceto P, Lai C, De Crescenzo F, et al. Cognitive decline after carotid endarterectomy: Systematic review and meta-analysis. *Eur J Anaesthesiol*. 2020;37(11):1066-1074.
13. Kordzadeh A, Abbassi OA, Prionidis I, Shawish E. The Role of Carotid Stump Pressure in Carotid Endarterectomy: A Systematic Review and Meta-Analysis. *Ann Vasc Dis*. 2020;13(1):28-37.
14. Cheng SF, Richards T, Gregson J, et al. Long Term Restenosis Rate After Carotid Endarterectomy: Comparison of Three Surgical Techniques and Intra-Operative Shunt Use. *Eur J Vasc Endovasc Surg*. 2021;62(4):513-521.
15. Chisci E, Lazzeri E, Masciello F, et al. "Timing to Carotid Endarterectomy Affects Early and Long Term Outcomes Of Symptomatic Carotid Stenosis.". *Ann Vasc Surg*. 2021.
16. Kumar V, Ramachandran S, Sylaja PN, Pitchai S. Conventional Carotid Endarterectomy with Shunt versus Eversion Carotid Endarterectomy without Shunt does the Technique Influence the Outcome in Symptomatic Critical Carotid Stenosis. *Asian J Neurosurg*. 2021;16(2):321-325.
17. Li Q, Liu B, Zhao Y, et al. Echolucent carotid plaque is associated with restenosis after carotid endarterectomy. *J Neurosurg*. 2020;134(3):1203-1209.
18. Rychen J, Madarasz A, Murek M, et al. Management of postoperative internal carotid artery intimal flap after carotid endarterectomy: a cohort study and systematic review. *J Neurosurg*. 2021:1-8.

19. Squizzato F, Xodo A, Tagliavero J, et al. Early outcomes of routine delayed shunting in carotid endarterectomy for symptomatic patients. *J Cardiovasc Surg (Torino)*. 2021;62(6):573-581.
20. Dakour-Aridi H, Gaber MG, Khalid M, Patterson R, Malas MB. Examination of the interaction between method of anesthesia and shunting with carotid endarterectomy. *J Vasc Surg*. 2020;71(6):1964-1971.
21. Chung BH, Heo SH, Park YJ, Kim YW, Woo SY, Kim DI. Comparative Analysis Using Propensity Score Matching Analysis: Primary Closure versus Patch Angioplasty During Carotid Endarterectomy. *Ann Vasc Surg*. 2020;62:166-172.
22. Chongruksut W, Vaniyapong T, Rerkasem K. Routine or selective carotid artery shunting for carotid endarterectomy (and different methods of monitoring in selective shunting). *The Cochrane database of systematic reviews*. 2014;2014(6):CD000190.
23. Wilkinson JM, Rochester JR, Sivaguru A, Cameron IC, Fisher R, Beard JD. Middle cerebral artery blood velocity, embolisation, and neurological outcome during carotid endarterectomy: a prospective comparison of the Javid and the Pruitt-Inahara shunts. *Eur J Vasc Endovasc Surg*. 1997;14(5):399-402.
24. Gaunt ME. Transcranial Doppler: preventing stroke during carotid endarterectomy. *Ann R Coll Surg Engl*. 1998;80(6):377-387.
25. Hayes PD, Vainas T, Hartley S, et al. The Pruitt-Inahara shunt maintains mean middle cerebral artery velocities within 10% of preoperative values during carotid endarterectomy. *J Vasc Surg*. 2000;32(2):299-306.
26. Kim TY, Choi JB, Kim KH, Kim MH, Shin BS, Park HK. Routine Shunting is Safe and Reliable for Cerebral Perfusion during Carotid Endarterectomy in Symptomatic Carotid Stenosis. *Korean J Thorac Cardiovasc Surg*. 2012;45(2):95-100.
27. Andrási TB, Kindler C, Dorner E, Strauch J. Transverse small skin incision for carotid endarterectomy. *Ann Vasc Surg*. 2015;29(3):447-456.
28. Bennett KM, Scarborough JE, Cox MW, Shortell CK. The impact of intraoperative shunting on early neurologic outcomes after carotid endarterectomy. *J Vasc Surg*. 2015;61(1):96-102.
29. Katano H, Yamada K. Comparison of internal shunts during carotid endarterectomy under routine shunting policy. *Neurologia medico-chirurgica*. 2014;54(10):806-811.
30. Wiske C, Arhuidese I, Malas M, Patterson R. Comparing the efficacy of shunting approaches and cerebral monitoring during carotid endarterectomy using a national database. *J Vasc Surg*. 2018;68(2):416-425.
31. Yuksel V, Canbaz S, Ege T, Sunar H. Is it necessary to use an intraluminal shunt in symptomatic patients with contralateral carotid artery stenosis? *Acta chirurgica Belgica*. 2014;114(3):179-182.
32. Tendera M, Aboyans V, Bartelink ML, et al. ESC Guidelines on the diagnosis and treatment of peripheral artery diseases: Document covering atherosclerotic disease of extracranial carotid and vertebral, mesenteric, renal, upper and lower extremity arteries: the Task Force on the Diagnosis and Treatment of Peripheral Artery Diseases of the European Society of Cardiology (ESC). *Eur Heart J*. 2011;32(22):2851-2906.
33. Naylor AR, Ricco JB, de Borst GJ, et al. Editor's Choice - Management of Atherosclerotic Carotid and Vertebral Artery Disease: 2017 Clinical Practice Guidelines of the European Society for Vascular Surgery (ESVS). *Eur J Vasc Endovasc Surg*. 2018;55(1):3-81.
34. Ricotta JJ, Aburahma A, Ascher E, Eskandari M, Faries P, Lal BK. Updated Society for Vascular Surgery guidelines for management of extracranial carotid disease. *J Vasc Surg*. 2011;54(3):e1-31.
35. Brott TG, Halperin JL, Abbara S, et al. 2011 ASA/ACCF/AHA/AANN/AANS/ACR/ASNR/CNS/SAIP/SCAI/SIR/SNIS/SVM/SVS guideline on the management of patients with extracranial carotid and vertebral artery disease. *Stroke*. 2011;42(8):e464-540.
36. Eckstein HH, Kühnl A, Dörfler A, Kopp IB, Lawall H, Ringleb PA. The diagnosis, treatment and follow-up of extracranial carotid stenosis. *Dtsch Arztebl Int*. 2013;110(27-28):468-476.
37. AbuRahma AF, Avgerinos ED, Chang RW, et al. The Society for Vascular Surgery implementation document for management of extracranial cerebrovascular disease. *Journal of vascular surgery*. 2022;75(1s):26s-98s.
38. Wafa HA, Wolfe CDA, Emmett E, Roth GA, Johnson CO, Wang Y. Burden of Stroke in Europe: Thirty-Year Projections of Incidence, Prevalence, Deaths, and Disability-Adjusted Life Years. *Stroke*. 2020;51(8):2418-2427.
39. Aboyans V, Ricco JB, Bartelink MEL, et al. 2017 ESC Guidelines on the Diagnosis and Treatment of Peripheral Arterial Diseases, in collaboration with the European Society for Vascular Surgery (ESVS): Document covering atherosclerotic disease of extracranial carotid and vertebral, mesenteric, renal, upper and lower extremity arteries Endorsed by: the European Stroke Organization (ESO) The Task Force for the Diagnosis and Treatment of Peripheral Arterial Diseases of the European Society of Cardiology (ESC) and of the European Society for Vascular Surgery (ESVS). *Eur Heart J*. 2018;39(9):763-816.
40. Setacci C, Argenterì A, Cremonesi A, et al. Guidelines on the diagnosis and treatment of extracranial carotid artery stenosis from the Italian Society for Vascular and Endovascular Surgery. *J Cardiovasc Surg (Torino)*. 2014;55(1):119-131.
41. Demarin V, Lovrencić-Huzjan A, Basić S, et al. Recommendations for the management of patients with carotid stenosis. *Acta Clin Croat*. 2010;49(1):101-118.
42. Venermo M, Wang G, Sedrakyan A, et al. Editor's Choice - Carotid Stenosis Treatment: Variation in International Practice Patterns. *Eur J Vasc Endovasc Surg*. 2017;53(4):511-519.

43. Clouse WD, Boitano LT, Ergul EA, et al. Contralateral Occlusion and Concomitant Procedures Drive Risk of Non-ipsilateral Stroke After Carotid Endarterectomy. *Eur J Vasc Endovasc Surg.* 2019;57(5):619-625.
44. AbuRahma AF, Aygerinos ED, Chang RW, et al. Society for Vascular Surgery clinical practice guidelines for management of extracranial cerebrovascular disease. *Journal of Vascular Surgery.* 2022;75(1):4S-22S.
45. Muller MD, Lyrer PA, Brown MM, Bonati LH. Carotid Artery Stenting Versus Endarterectomy for Treatment of Carotid Artery Stenosis. *Stroke.* 2021;52(1):e3-e5.
46. Lee EJ, Cho YP, Lee SH, et al. Hemodynamic Tandem Intracranial Lesions on Magnetic Resonance Angiography in Patients Undergoing Carotid Endarterectomy. *Journal of the American Heart Association.* 2016;5(10).
47. Levin SR, Farber A, Goodney PP, et al. Shunt intention during carotid endarterectomy in the early symptomatic period and perioperative stroke risk. *J Vasc Surg.* 2020;72(4):1385-1394 e1382.
48. Asensio JA, Kessler JJ, 2nd, Kotaru TR, Kalamchi LD, Miljkovic SS, Dabestani PJ. Penetrating Internal and Common Carotid Artery Injuries Shunts versus no shunts during repair effect on neurological outcomes. *Injury.* 2021;52(2):266-273.
49. Yang SS, Kim YW, Kim DI, et al. Impact of contralateral carotid or vertebral artery occlusion in patients undergoing carotid endarterectomy or carotid artery stenting. *J Vasc Surg.* 2014;59(3):749-755.
50. Kondov S, Beyersdorf F, Schöllhorn J, et al. Outcome of Near-Infrared Spectroscopy-Guided Selective Shunting During Carotid Endarterectomy in General Anesthesia. *Ann Vasc Surg.* 2019;61:170-177.
51. Inčiūra D, Antuševas A, Aladaitis A, Gimžauskaitė A, Velička L, Kavaliauskienė Ž. Near-infrared spectroscopy as a predictor of cerebral ischaemia during carotid endarterectomy in awake patients. *Vascular.* 2020;28(3):301-308.
52. Hicks KA, Stockbridge NL, Targum SL, Temple RJ. Bleeding Academic Research Consortium consensus report: the food and drug administration perspective. In. Vol 123: Am Heart Assoc; 2011:2664-2665.
53. Antuševas A, Aladaitis A, Velička L, et al. Outcomes of simultaneous carotid endarterectomy and coronary artery bypass grafting: A single centre experience. *Vascular.* 2023;31(5):914-921.
54. Grillo VTRS, Jaldin RG, Bertanha M, Sobreira ML, Soares CSP, de Camargo PAB. What to do when advanced thyroid cancer invades the carotid artery? Therapeutic challenge. *Jornal Vascular Brasileiro.* 2022;21.
55. Lee J, Lee S, Kim SW, Chang JW. Selective Shunting Based on Dual Monitoring with Electroencephalography and Stump Pressure for Carotid Endarterectomy. *Vasc Specialist Int.* 2018;34(3):72-76.

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