Reddick[®]-Saye Screw Retractor Kit

Instructions for Use – English

Reddick[®]-Saye Wundhaken-Kit

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Kit divaricatore a vite Reddick[®]-Saye

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Kit retractor de rosca Reddick[®]-Saye

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Kit Retractor Helicoidal Reddick[®]-Saye

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STERILE EO

Description

The Reddick®-Saye Screw Retractor Kit is specifically designed for use during laparoscopic cholecystectomy procedures. The retractor is 30 cm in length and 5 mm in diameter with allows insertion through standard trocars. The cork screw tip grasps the diseased gallbladder for manipulation. The Reddick-Saye needle set is designed to puncture the gallbladder and allow easy removal of the bile.

Instructions

The Reddick-Saye Retractor is utilized to grasp the fundus of the acute or chronically inflamed gallbladder. The gallbladder wall is thickened in these diseased processes, making it difficult or impossible to hold the fundus with standard laparoscopic forceps. The Reddick-Saye needle and sheath are inserted percutaneously into the fundus of the gallbladder. The needle is removed leaving the sheath inside the gallbladder. The bile is drained with syringe aspiration. The sheath is removed from the abdomen. The Reddick-Saye Retractor is introduced through the lateral laparoscopic portal and aimed toward the fundus of the gallbladder. The tip of the Reddick-Saye Retractor is introduced into the thickened gallbladder wall and the Reddick-Saye Retractor turned clockwise two revolutions. After adequate anchoring into the wall, the push rod is advanced superiorly to push the gallbladder and liver toward the diaphragm.

A second Reddick-Saye Retractor can be placed in the Hartmann's pouch if it is also thickened. The Reddick-Saye Retractor is introduced through the midclavicular port. The Reddick-Saye Retractor is turned clockwise two revolutions. Once anchoring has occurred, the push rod is pulled inferiorly and laterally to expose the cystic duct and artery.

At the conclusion of the procedure, both Reddick-Saye Retractors are turned two revolutions counter clockwise and removed through the laparoscopic ports. The gallbladder is removed through the umbilical portal with either a standard grasper or a Reddick-Saye Retractor.

WARNING

DO NOT INSERT THE REDDICK-SAYE RETRACTOR MORE THAN TWO REVOLUTIONS AS THE SHARP TIP MAY DAMAGE SURROUNDING TISSUE.

WARNING

DRAIN THE BILE BEFORE INSERTING THE REDDICK-SAYE RETRACTOR TO REMOVE THE GALLBLADDER.

Cautions

1. The Reddick-Saye Retractor is for temporary use only
2. The Reddick-Saye Retractor is for single use only.
3. The Reddick-Saye Retractor is supplied sterile and non-pyrogenic in a peel-open package.
4. Inspect each package and do not use if there is any evidence that the package has been opened or damaged.
5. Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

Notices: Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and Reddick are registered trademarks of LeMaitre Vascular, Inc. U.S. Patent #5,437,266 and #5,573,496.
©2013 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH