## **Martin™ Dissector**

Instructions for Use - English

## **Martin™ Dissector**

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## **Dissector Martin™**

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## **Dissector Martin™**

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**Martin™ Dissector**

**Martin™ Dissector**  
(Model Numbers 4200-10, e4200-10)  
Remote Endarterectomy Device  
Instructions for Use - English

**STERILE EO** 

**Device Description**

The Martin™ Dissector consists of a shaft with an angled spatula at the distal end. The shaft and angled spatula are constructed of stainless steel. The inner lumen of the shaft extends from the handle to the palm of the spatula.

The Martin Dissector is available in one size. The useable length of the spatula/shaft is approximately 40 cm and the hollow lumen can accommodate up to a 0.018" diameter guidewire.

**Indications And Usage**

The Martin Dissector is intended for remote intraluminal stripping of peripheral blood vessels during vascular reconstruction and to provide a conduit for guidewire placement.

General indications for peripheral vessel reconstruction include disabling claudication, ischemic pain at rest, ischemic skin lesions or gangrene.

Indications for endarterectomy specifically as a procedure for vessel reconstruction include localized lesions and heavy circumferential plaque deposits.

**Contraindications**

Contraindications include, but are not limited to:

- Any evidence of aneurysmal change or severe calcification may contraindicate endarterectomy since continued aneurysmal degeneration of the adventitia may be expected.
- In some cases, the vessel wall layers may be quite coherent and the plaque segment may be difficult to separate from the vessel wall. When the vessel is in this state of disease, it is often advisable to abandon endarterectomy and consider an alternative method of vascular reconstruction.
- Use of this device is contraindicated in the coronary, internal carotid arteries and other vessels of the central circulatory system.

**Warnings And Precautions**

The Martin Dissector may only be used by physicians trained in vascular surgery, specifically endarterectomy procedures. The techniques described herein are offered to clinicians as a guide.

**Martin Dissector - Handling Precautions**

Always inspect the sterile packaging before use. If the packaging is damaged, do not use the product and contact your LeMaitre Vascular representative. Do not re-sterilize the Martin Dissector. The Martin Dissector must be used prior to the "Use By" date printed on the package.

**Martin Dissector - Use Precautions**

Always use fluoroscopy when advancing or manipulating endovascular devices. Do not attempt to forcibly advance the Martin Dissector.

Use of the Martin Dissector may lead to vessel perforation requiring additional interventional procedures (e.g., bypass graft surgery of the artery).

**Potential Adverse Events**

Adverse events potentially associated with the use of the Martin Dissector (e.g., all complications associated with PTA, remote endarterectomy or atherectomy procedures) include but are not limited to the following:

- Incisional/puncture site inflammation, hematoma or bleeding
- Wound infection (incisional)
- Vessel injury, dissection, perforation, rupture, or pseudoaneurysm
- Vessel spasm
- Emboli, thrombosis or acute occlusion
- Fever, drug reactions or allergic reactions to contrast medium

- Conversion to bypass
- Limb loss due to distal embolization
- Renal failure
- Death

## **Contents**

One (1) Sterile Martin Dissector Device One (1) Instructions for Use

STORAGE: Store in a cool, dry place.

INSPECTION: Carefully inspect the sterile package for damage or defects before opening. Do not use product after the "Use By" date on the package. If the integrity of the sterile package has been compromised or the package or product is defective, do not use the product and contact LeMaitre Vascular for return information.

## **Recommended Material**

Prepare the following material using sterile technique:

- Up to a 0.018" guidewire of appropriate length (recommended guidewire: TSCF-18-145-3 by Cook Inc.)

## **Directions For Use**

Specific patient preparation and procedural techniques are determined by the physician. One suggested technique is as follows:

1. Patient preparation and administration of anticoagulation/antiplatelet regimen is left to the discretion of the physician.
2. Use sterile surgical precautions and standards.
3. Using appropriate anesthesia, gain access at the appropriate site through an arteriotomy to expose the atheroma.
4. Make a slanting, bevel-type incision of the proximal atheroma to start a cleavage plane between the atheroma and surrounding media.

*Note: The following procedure steps should be performed under fluoroscopic guidance.*

5. Insert the Martin Dissector into the cleavage plane of the vessel. Using a gentle pushing motion around the atheroma core, advance the device to a point distal to the occluded segment.
6. Under fluoroscopic guidance, position the Martin Dissector at the end of the dissection plane.
7. An 0.018" guidewire maybe inserted into the Martin Dissector lumen and advance until the distal tip of the wire is sufficiently anchored into the true vessel lumen. Anchoring of the guidewire will allow retraction of the Martin Dissector while maintaining guidewire position.

*Note: Do not use polymer covered hydrophilic guidewires.*

*Optimum guidewire performance is achieved with the use of flex wires or wires with very long flexible tips.*

8. Remove the Martin Dissector from the patient while maintaining guidewire position.
9. Following completion of the remote endarterectomy procedure, maintain guidewire position to facilitate additional interventional procedures as required.
10. Close arteriotomy using standard procedure.

## **Resterilization/Re-use**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

### **Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed accorded to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

#### **Cleaning:**

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

#### **Packaging:**

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803

### **Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this

device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre is a registered trademark, and Martin is a trademark,
©2020 LeMaitre Vascular, Inc. All rights reserved. Printed in
LeMaitre Vascular GmbH
of LeMaitre Vascular, Inc.