## **LeMills Valvulotome**

Instructions for Use - English

## **LeMills Valvulotom**

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## **Valvulotome LeMills**

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## **Valvulotomo LeMills**

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## **Valvulotomo LeMills**

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## **Valvulótomo LeMills**

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**STERILE** | **EO** | **Rx only**

### Description

The LeMills family of Valvulotomes are devices that cut venous valves during vascular procedures. The devices consist of stainless steel blades and plastic handles.

### Indication for Use

The LeMills Valvulotomes are intended to cut venous valves.

### Contraindications

1. endarterectomy procedures,
2. thrombolysis procedures,
3. vein stripping procedures, or
4. embolectomy procedures.

### Warnings

1. Do not use if the package or device is damaged.
2. Do not use if inner packaging is opened outside a sterile environment.
3. The LeMills Valvulotomes are to be used by a qualified physician.
4. Do not pass the LeMills Valvulotome through a vessel that has undergone synthetic grafting or contains implants.
5. Do not use device for valvulotomy unless the target vein is fully distended by arterial blood flow (due to anastomosis of artery to vein) or by saline injection.

### Precautions


1. United States Federal and other law restricts this device to sale on or by the order of a physician.
2. Do not use devices past the expiration date printed on the labeling.
3. These devices are single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/Repackaging" section of this document for further details). The device is NOT autoclavable.
4. Keep any dissected portion of the vein straight and avoid twisting the vein during valvulotomy.
5. Use caution when tracking the LeMills Valvulotomes in thrombophlebitic veins.

### Potential Complications


- Vessel wall perforation
- Intimal disruption
- Entrapment of cutting blades in branches or anastomoses
- Inadequate valve disruption
- Vein Stricture
- Thrombus formation
- Post-procedure wound necrosis

### Specifications

Antegrade LeMills Valvulotome:

Retrograde LeMills Valvulotome:

Catalog Number	Description	Useable Length (cm)
1050-00, e1050-00	Retrograde LeMills Valvulotome	12.5
1050-01, e1050-01	Antegrade LeMills Valvulotome	12.5
1050-02, e1050-02	Antegrade & Retrograde LeMills Valvulotomes	12.5; 12.5

### To Open Package

1. Open box and remove sealed sterile unit.
2. Open sealed sterile unit and present contents to personnel in sterile environment.

### Pre-Use Check

Inspect blades for damage.

### Instructions for Use

1. Insert the device into the vein to the treatment area.
2. If using Retrograde LeMills, slowly retract the device to disrupt the valve sets. If using Antegrade LeMills, slowly advance the device to disrupt the valve sets.
3. Steps 1 and 2 may be repeated as necessary.
4. Confirm the free flow of blood through the vein.

### Storage

Store in a cool, dry place.

### Resterilization/Re-use

These devices are single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

### Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre is a registered trademark of LeMaitre Vascular, Inc.
©2014 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH