**Pruitt F3®-S Polyurethane Carotid Shunt**

Instructions for Use - English

**Pruitt F3®-S Polyurethan-Karotisshunt**

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#### Introduction

The Pruitt F3<sup>®</sup>-S Polyurethane Carotid Shunt is designed to serve as an artificial passage connecting two blood vessels, allowing blood flow from one vessel to another. This is accomplished by using a clear, plastic, sterile conduit that is held in place by a stabilization technique on both ends of the conduit.

#### Product Description

The Pruitt F3<sup>®</sup>-S Polyurethane Carotid Shunt (the Shunt) is a multi-lumen device with balloons at both the distal (internal carotid) and proximal (common carotid) ends of the Shunt. The balloons, when inflated independently, act as a stabilization mechanism to maintain the position of the Shunt when it is placed within the common and internal carotid arteries. The Pruitt F3<sup>®</sup>-S Polyurethane Carotid Shunt has features to aid the user during Shunt insertion and balloon inflation. The inflation path of the proximal (common carotid) balloon is color-coded. Sterile saline is injected from the blue stopcock, through the blue lumen and into the blue common carotid balloon. Depth markings on the shunt body are for reference during insertion.

#### Indication

1. Carotid Shunts are indicated for use in carotid endarterectomy as a temporary conduit to allow for blood flow between the common and internal carotid arteries.

#### Contraindications

1. The Shunt is a temporary device and should not be implanted.
2. The Shunt is not indicated for use in embolectomy, thrombectomy, or vessel dilation.

#### Warnings

1. Do not reuse. Do not resterilize. For single use only.
2. Do not use air or gas to inflate the balloons. Inflate the balloons with sterile saline.
3. Do not inflate the internal carotid balloon to any greater volume than is necessary to obstruct blood flow for the internal carotid artery. DO NOT EXCEED the recommended maximum balloon liquid capacity (see Specifications).
4. Exercise caution when encountering extremely diseased vessels. Arterial rupture or balloon failure due to sharp calcified plaque may occur.
5. Deflate the balloons prior to Shunt removal. Avoid using excessive force to push or pull the Shunt against resistance.

#### Precautions

1. Inspect the product and package prior to use and do not use if there is any evidence that the package or the Shunt has been damaged.
2. The Shunt should be used only by qualified physicians thoroughly familiar with cardiovascular surgical procedures involving the carotid artery. Described procedures are provided for informational purposes only. Each physician must determine the appropriate use of this device for each patient.
3. Pretest the Shunt according to the pretest procedure prior to patient use to ensure the lumen is free of obstructions and the balloons are functional.
4. Aspirate the balloons prior to inflation.
5. Place clear balloon into internal carotid artery and blue balloon in common carotid artery.
6. If the Shunt is not properly maintained in position through balloon stabilization, it may migrate within the internal carotid artery, potentially damaging the intima.
7. Avoid extended or excessive exposure to fluorescent light, heat, sunlight, or chemical fumes to reduce balloon degradation. Excessive handling during insertion, and/or plaque and other deposits within the blood vessel, may damage the balloon and increase the possibility of balloon rupture.
8. Do not grasp the balloon with instruments at any time to avoid damage to the balloon.
9. Make secure connections between the syringe and the hub to avoid introduction of air.
10. Dispose of used product in accordance with established hospital protocols for biohazards.
11. Do not reinsert the Shunt once it has been removed, it may become occluded from stagnant blood left in the Shunt.

#### Adverse Events

As with all cardiovascular procedures involving the carotid arteries, complications may occur during or following carotid endarterectomy. These may include, but are not limited to:

- stroke
- transient ischemic attack
- neurologic complications

- embolization of blood clots, arteriosclerotic plaque, or air
- hypertension or hypotension
- infection
- intimal disruption
- arterial dissection
- vessel perforation and rupture
- hemorrhage
- arterial thrombosis
- aneurysms
- arterial spasm

**How Supplied**

The Shunt is supplied sterile. The sterility of the package is assured as long as it is unopened and undamaged.

**Procedure**

**IMPORTANT:** A variety of surgical techniques may be used when using the Shunts; therefore, the surgeon is best advised to use the method which his/her own practice and training dictate to be best for the patient. Specific surgical techniques are left to the discretion of the surgeon.

**Balloon Pretest Procedure (Perform Before Patient Use)**

1. Inflate both balloons up to the maximum recommended volumes with sterile saline and inspect for leaks. If there is any evidence of leaks around the balloons or if either balloon will not remain inflated, do not use the product.
- NOTE: The common carotid balloon is designed to inflate partially to minimize pressure on the common carotid artery while maintaining position.*
2. Before patient use, aspirate the balloons completely prior to inflation of the balloons.

**T-Port Pretest (Perform Before Patient Use)**


1. Place a gloved finger over the opening near the large blue balloon end and inject sterile saline through the red T-Port stopcock. Fluid should flow through the opening near the small clear balloon end.
2. Place a gloved finger over the small clear balloon end and inject sterile saline through the red T-Port stopcock. Fluid should flow through the opening near the large blue balloon end.
3. Do not use the Shunt if fluid does not flow through both openings.

**Procedure for T-Port Models 2016-10, 2016-11, 2014-10, 2014-12, 2015-10**

1. Expose the carotid artery and perform the arteriotomy in the usual manner.
2. Place a clamp on the proximal (common carotid) end of the Shunt between the blue balloon and the T-Port arm.
3. Place the distal end (clear balloon) of the Shunt into the internal carotid artery.
4. Attach a 3 ml syringe to the white stopcock and SLOWLY inflate the clear balloon with up to 0.25 ml of sterile saline. (Figure A)
5. As inflation progresses, carefully observe back-bleeding from the internal carotid artery around the Shunt. The back-bleeding will diminish as the balloon expands. When the balloon is inflated sufficiently to occlude the artery, back-bleeding around the Shunt will stop, there will be a feeling of slight resistance to further inflation. This is the end-point: STOP INFLATION IMMEDIATELY AT THIS POINT. (Figure B)
6. Close the white stopcock.
7. Open the T-Port stopcock and allow blood to back-bleed through the T-Port of the Shunt observing for air bubbles and/or atheromatous debris.
8. When no debris or bubbles are noted, close the T-Port stopcock. Move the clamp to the distal (internal carotid) end of the Shunt between the white balloon and the T-port arm.
9. Place the proximal (blue balloon) end of the Shunt into the common carotid artery.
10. Attach a 3 ml syringe to the blue stopcock and slowly inflate the blue common carotid artery balloon with up to 1.5 ml of sterile saline. Close the blue stopcock.
11. Open the T-Port stopcock and allow blood to flow through the T-Port of the Shunt observing for air bubbles and/or atheromatous debris.
12. When no debris or bubbles are noted, close the T-Port stopcock and remove the clamp on the distal (internal carotid) end of the Shunt. Proceed with the procedure.
13. When the endarterectomy is completed, deflate the balloons, remove the Shunt and close the arteriotomy in the usual manner. (Figure C)

**Procedure for Models 2014-11 and 2014-13**

1. Expose the carotid artery and perform the arteriotomy in the usual manner.
2. Place the distal (clear balloon) end of the Shunt into the internal carotid artery.

Appearance of balloon in Internal Carotid Artery


Figure A

Figure B

Figure C

3. Attach the 3 ml syringe to the white stopcock and SLOWLY inflate the clear balloon with up to 0.25 ml of sterile saline. (Figure A)
4. As inflation progresses, carefully observe back-bleeding from the internal carotid artery around the Shunt. The back-bleeding will diminish as the balloon expands. This is the end-point: STOP INFLATION IMMEDIATELY AT THIS POINT. (Figure B)
5. Close the white stopcock.
6. Place the proximal (blue balloon) end of the Shunt into the common carotid artery.
7. Attach a 3 ml syringe to the blue stopcock and slowly inflate the blue common carotid artery balloon with up to 1.5 ml of sterile saline, close the blue stopcock.
8. Clamp across the Shunt lumen and remove the clamp on the common carotid artery. Slowly remove the clamp that is across the Shunt and observe for air bubbles and/or atheromatous debris. If no bubbles or debris are seen, fully remove the clamp. Proceed with the procedure.
9. When the endarterectomy is completed, deflate the balloons, remove the Shunt and close the arteriotomy in the usual manner. (Figure C)

**Specifications**

Model	Description	Usable Length	Diameter	Common Inflation Lumen Markings
2016-10	Pruitt F3-S Polyurethane Carotid Shunt with T-port (Outlying)	31 cm	10 French (3.3 mm)	Blue Lumen
2016-11	Pruitt F3-S Polyurethane Carotid Shunt with T-port (Inlying)	15 cm	10 French (3.3 mm)	Blue Lumen
2014-10	Pruitt F3-S Polyurethane Carotid Shunt with T-Port (Outlying)	31 cm	9 French (3.0 mm)	Blue Lumen
2014-11	Pruitt F3-S Polyurethane Carotid Shunt (Outlying)	31 cm	9 French (3.0 mm)	Blue Lumen
2014-12	Pruitt F3-S Polyurethane Carotid Shunt with T-Port (Inlying)	15 cm	9 French (3.0 mm)	Blue Lumen
2014-13	Pruitt F3-S Polyurethane Carotid Shunt (Inlying)	15 cm	9 French (3.0 mm)	Blue Lumen
2015-10	Pruitt F3-S Polyurethane Carotid Shunt with T-Port (Outlying)	31 cm	8 French (2.67 mm)	Blue Lumen

	Stopcock Color	Balloon Maximum Liquid Capacity
Common Carotid Balloon	Blue	1.5 ml
Internal Carotid Balloon	White	.25 ml
T-Port	Red	N/A

		Internal Balloon Volume (mL)				
		0.05 mL	0.10 mL	0.15 mL	0.20 mL	0.25 mL
Average Balloon Diameter (mm)	8F	4.7	5.8	6.7	7.4	8.0
	9F	5.5	6.6	7.3	7.9	8.4
	10F	5.3	6.1	6.9	7.6	8.3

		Common Balloon Volume (mL)					
		0.25 mL	0.50 mL	0.75 mL	1.00 mL	1.25 mL	1.50 mL
Average Balloon Diameter (mm)	8F	6.7	9.3	10.9	12.1	12.9	13.7
	9F	7.3	9.6	11.1	12.3	13.2	14.0
	10F	6.6	9.3	10.9	12.2	13.3	14.0

**Storage Shelf Life**

The shelf life is indicated by the USE BY date on the package label. The USE BY date printed on each label is NOT a sterility date. The USE BY date is based on the normal life expectancy of the Shunt when properly stored. The use of the Shunt beyond the expiration date is not recommended because of potential bond deterioration. LeMaitre Vascular, Inc. does not make provisions for replacing or reprocessing expired product.

Proper storage procedures must be practiced to achieve optimum shelf life. The product should be stored in a cool dark area. Proper stock rotation should be practiced.

**Restoration/Re-use**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

**Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

**Pruitt F3®-S Polyurethan-Karotisshunt (außenliegend)** (Modell-Nr. 2016-10, 2014-10, 2014-11, 2015-10)  
**Pruitt F3®-S Polyurethan-Karotisshunt (innenliegend)** (Modell-Nr. 2016-11, 2014-12, 2014-13)  
Gebrauchsanweisung – Deutsch

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LeMaitre, Pruitt, and Pruitt F3 are registered trademarks of LeMaitre Vascular, Inc.
©2018 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH